## Adenosina Deaminase - ADA

Adenosine Deaminase / Adenosina Deaminasa  
Ref. 11.063.00

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#### ENGLISH

#### INTENDED USE

Kit intended to determine adenosine deaminase (ADA) in serum, plasma, pleuroperitoneal fluid and cerebrospinal fluid samples. Diagnostic use only.




#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- It is recommended to calibrate the product every two weeks.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Enzymatic – Peroxidase

Adenosine deaminase (ADA) promotes the enzymatic deamination of adenosine to inosine, which is converted to hypoxanthine by purine nucleoside phosphorylase (PNP). Xanthine oxidase (XOD) catalyzes the conversion of hypoxanthine to uric acid and hydrogenous peroxide, which reacts with TOOS and 4-aminoantipyrine in the presence of peroxidase (POD), to generate the quinoneimine dye. The last one is spectrophotometrically measured at 546 nm. Color intensity is proportional to the concentration of adenosine deaminase in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum, plasma (EDTA or heparin), pleuroperitoneal fluid and cerebrospinal fluid.

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Serum and Plasma	2 to 8 °C	2 weeks
	-20 °C	6 months
Cerebrospinal Fluid	2 to 8 °C	1 week
	-20 °C	1 month
Pleuroperitoneal Fluid	-20 °C	6 days

#### PRODUCT DESCRIPTION

**R 1** Tris-HCl buffer 25 mmol/L pH 8,0; peroxidase 100 U/L; 4-aminoantipyrine 4 mmol/L; PNP 300 U/L; XOD 400 U/L

**R 2** Tris-HCl buffer 25mmol/L pH 4,0; adenosine 11,0 mmol/L; TOOS 4,0 mmol/L; preservative.

**CAL** Hepes buffer 50 mmol/L; adenosine deaminase 50 U/L. Traceable to reference material BCR-647.

**CONTROL 1** Hepes buffer 50 mmol/L; adenosine deaminase  $\cong$  30 U/L

**CONTROL 2** Hepes buffer 50 mmol/L; adenosine deaminase  $\cong$  130 U/L

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratorial quality control, it is recommended the use of the control below:

Calibrador – ADA CAL	XX.XXX.XX
Level 1 Control – ADA CONTROL 1	XX.XXX.XX
Level 2 Control – ADA CONTROL 2	XX.XXX.XX

**REF**

#### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 546 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

##### A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

Reaction Type	Kinetic / Positive
Wavelength	546 nm
Temperature	37 °C
Sample volume*	5 µL
R1 Volume*	180 µL
Sample + R1 Incubation Time	300 seconds
R2 Volume*	90 µL
R1 + Sample + R2 Incubation Time	210 seconds
Measure the absorbance continuously for 90 seconds	
Calibration	2 points

\* Volumes can be modified as long as the established proportion is maintained.

##### B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

##### C) INTERPRETATION

Adenosine deaminase (ADA) is widely distributed in several tissues of the human body, being found mainly in T lymphocytes. Increased activity of serum adenosine deaminase is observed in acute hepatitis, alcoholic liver fibrosis, chronic active hepatitis, cirrhosis, viral hepatitis and it can be observed in tuberculosis exudate. Combined with ALT or γ-GT. The detection of serum adenosine deaminase activity in patients is valuable for the diagnosis of hepatic disease and its determination in body fluid samples can help to differentiate the diagnosis of tuberculosis.

#### PERFORMANCE CHARACTERISTICS

Operating range
3,0 a 200,0 U/L

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0,88 U/L	1,17 U/L

Analytical Specificity		
Ascorbic Acid	Hemoglobin	Triglycerides
50 mg/dL	800 mg/dL	1000 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y = 1,0224x + 0,0066
Correlation Coefficient (R)	0,9994

The product's accuracy was evaluated with the reference material BCR-647, at the decision level of 30,6 U/L, with a total systematic error of 1,53%.

#### Precision:

Within-run precision was determined with 20 replicates in an analytical run; between-run precision was determined with 5 replicates for 5 days.

Samples (U/L)	Within-Run Precision		Between-Run Precision	
	SD (U/L)	CV (%)	SD (U/L)	CV (%)
30.0	0.170	0.59	0.030	1.70
140.0	0.286	0.20	0.110	0.51

CV: Coefficient of variation; SD: Standard deviation.

#### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

#### REFERENCE RANGES

Serum and Plasma	0 - 15 U/L
Pleuroperitoneal Fluid	0 - 30 U/L
Cerebrospinal Fluid	0 - 9 U/L

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

## WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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- UNGERER, J. P. et al. Serum adenosine deaminase: isoenzymes and diagnostic application. *Clinical chemistry*, v. 38, n. 7, p. 1322-1326, 1992.
- SHI, P.; LI, W.; LI, Y. The clinical significance of serum adenosine deaminase determination in chronic hepatitis B patients [J]. *Laboratory Medicine*, v. 6, 2009.
- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AAC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

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