## Amônia - AMN

Ammonia / Amônia  
Ref. 12.015.00

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### ENGLISH

#### INTENDED USE

Kit intended to determine ammonia in plasma samples. Diagnostic use only.


#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- It is recommended to calibrate the product every 3 days.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8 °C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** Enzymatic

Ammonia from the sample, alongside NADPH and  $\alpha$ -ketoglutaric acid, is catalyzed by glutamate dehydrogenase (GLDH) to form glutamic acid and NADP<sup>+</sup>. The concentration of ammonia in the sample is determined by the consumption rate of NADPH, which is spectrophotometrically measured at 340 nm.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** plasma EDTA.

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

**Preservation:**

	Temperature	Stability Period
Plasma	4 to 8 °C	4 hours
	-20 °C	3 weeks

### PRODUCT DESCRIPTION

**R 1** NADPH > 0,1 mmol/L and  $\alpha$ -ketoglutaric acid > 5 mmol/L.

**R 2** Glutamate dehydrogenase > 0,5 kU/L

**CAL** Buffer 100 mmol/L and ammonium chloride in the concentration range of 200 – 400 mmol/L.

**CONTROL 1** Buffer 100 mmol/L and ammonium chloride in the concentration range of 50 – 300 mmol/L.

**CONTROL 2** Buffer 100 mmol/L and ammonium chloride in the concentration range of 250 – 500 mmol/L.

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrator – Ammonia CAL

Level 1 Control – Ammonia CONTROL 1

Level 2 Control – Ammonia CONTROL 2

<b>REF</b>	XX.XXX.XX
	XX.XXX.XX
	XX.XXX.XX

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 340 nm.
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

This product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnicaitda.ind.br](http://www.biotecnicaitda.ind.br)

Reaction Type	Endpoint / Negative
<b>Wavelength</b>	340 nm
<b>Temperature</b>	37 °C
<b>Sample volume*</b>	17 $\mu\text{L}$
<b>R1 Volume*</b>	150 $\mu\text{L}$
<b>Sample + R1 Incubation Time</b>	300 seconds
<b>R2 Volume*</b>	50 $\mu\text{L}$
<b>R1 + Sample + R2 Incubation Time</b>	300 seconds
<b>Calibration</b>	2 points

\* Volumes can be modified as long as the established proportion is maintained.

### B) CALCULATIONS

The equipment automatically calculates the analyte's concentration in each sample.

### C) INTERPRETATION

Under normal circumstances, ammonia is converted to urea in the liver. In cases of severe hepatic disease, the conversion of ammonia into urea does not happen; therefore, ammonia is not eliminated from the blood and there is an increase in its concentration. When in elevated concentrations is the blood, ammonia is neurotoxic, which may lead to hepatic encephalopathy (hepatic coma). The determination of plasmatic ammonia in adults is mainly used to monitor and treat hepatic coma. Besides, the determination of plasmatic ammonia can be used in the pediatric diagnosis of Reye's syndrome.

### PERFORMANCE CHARACTERISTICS

Operating range
10.65 to 1000.00 $\mu\text{mol/L}$

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
10.88 $\mu\text{mol/L}$	10.65 $\mu\text{mol/L}$

Analytical Specificity		
Ascorbic Acid	Bilirubin	Triglycerides
25 mg/dL	10 mg/dL	100 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 samples
Regression Equation	y = 0.9958x - 0.5057
Correlation Coefficient (R)	0.9987

The product's accuracy was evaluated with a sample spiked with ammonium chloride, at the decision levels of 109.1  $\mu\text{mol/L}$  and 218.2  $\mu\text{mol/L}$ , with a total systematic error of 2,5% and 1,2%, respectively.

### Precision:

Within-run precision was determined with 20 replicates in an analytical run; between-run precision was determined with 3 replicates, in two periods, for 20 days.

Samples ( $\mu\text{mol/L}$ )	Within-Run Precision		Between-Run Precision	
	SD ( $\mu\text{mol/L}$ )	CV (%)	SD ( $\mu\text{mol/L}$ )	CV (%)
120.00	2.18	1.78	4.67	3.89
435.00	2.10	0.48	6.15	1.41

CV: Coefficient of variation; SD: Standard deviation.

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

### REFERENCE RANGES

Plasma	0 - 47 $\mu\text{mol/L}$
--------	--------------------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotecnicaitda.com.br](mailto:sac@biotecnicaitda.com.br).

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- Department of Medical Administration, Ministry of Health of the People's Republic of China. National Clinical Laboratory Procedures [M]. 3<sup>rd</sup> edition, Nanjing: Southeast University Press, 2006, 460-462.

- YOUNG, D.S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. Clin. Chem. v.27 p.493-501, 1981.

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