COBRE

COPPER / COBRE
Ref. 12.013.00

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ENGLISH

INTENDED USE

Kit intended to determine copper ions in serum and heparin plasma samples. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- It can be observed a color alteration in R1 from transparent to brown, however it does not affect the reagent's reactional capacity as long as it is properly calibrated.
- A precipitate can be observed in R1 when stored under intense refrigeration. This does not affect the reagent's reactional capacity as long as it reaches room temperature before being used. Homogenize the product by gentle inversion before use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: di-Br-PAESA

Around 95% of body's copper ions is found bounded to ceruloplasmin, a protein responsible for the regulation of iron transport through the organism's membranes. In acidic pH, ions copper dissociate from ceruloplasmin and react with the chromogen di-Br-PAESA, forming a colored complex that can be measured at 578 nm (535-605 nm). Color intensity is proportional to the concentration of ions copper in sample.

SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

Sample Type: serum and heparin plasma.

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Serum and Plasma	4 to 8 °C	2 weeks
	-20 °C	1 year

PRODUCT DESCRIPTION

R 1	Acetate buffer < 500 mM; surfactants and stabilizers.
R 2	3,5 Di-Br-PAESA pH > 3,0
STD	Copper solution at concentration of 220 µg/dL and stabilizer. Traceable to reference materials NIST-SRM 3114 and NIST-SRM 928.

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the control below:

Pathological Control – Quantialt	REF	13.004.00
----------------------------------	------------	-----------

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 578 nm (535 – 605 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Pre-analytical recommendations:

- Let the reagents reach room temperature. Homogenize by gentle inversion before use.
- Wash the equipment before its use. Preferably, perform the Copper assay isolated from the routine.

1. Pipette in the test tubes:

	Blank	Sample Blank	Sample	STD
STD	---	---	---	60 µL
Sample	---	60 µL	60 µL	---
Purified water	60 µL	---	---	---
R1	900 µL	900 µL	900 µL	900 µL
Homogenize and incubate the test tubes for 3 minutes at 37°C				
R2	100 µL	---	100 µL	100 µL

2. Homogenize and incubate at 37 °C for 3 minutes.

3. Measure the Blank's, Sample Blank's, Sample's and Standard's absorbance at 578 nm (535-605 nm). The final reaction is stable for 30 minutes.

B) CALCULATIONS

Copper (µg/dL) = $\frac{\text{Sample's } \Delta A - \text{Blank} \times \text{STD Concentration} (\mu\text{g/dL})}{\text{STD's } \Delta A - \text{Blank}}$

In which:

ΔA = Sample's or STD's Absorbance – Sample or STD Blank's Absorbance

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration} (\mu\text{g/dL})}{\text{STD's } \Delta A - \text{Blank}}$

Copper (µg/dL) = $(\text{Sample's } \Delta A - \text{Blank}) \times \text{Calibration Factor}$

Automação: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaltda.ind.br

C) INTERPRETATION

Ions copper are cofactor to various enzymatic reactions, besides being part of iron transport and absorption metabolism. Reduced levels of copper are correlated with diminished serum proteins, which is observed in the following cases: malnutrition, incorrect absorption (celiac disease), loss of protein through feces or urine (nephrotic syndrome) or Wilson's disease. Increased levels of copper can be observed during pregnancy, chronic or acute infections, myocardial infarction, hyperthyroidism and hematological diseases.

PERFORMANCE CHARACTERISTICS

Operating range
15,3 a 640,0 µg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
9,60 µg/dL	9,80 µg/dL

Analytical Specificity			
Hemoglobin	Bilirubin	Triglycerides	Uric Acid
500 mg/dL	35 mg/dL	2000 mg/dL	40 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	$y = 1,0289x - 3,4441$
Correlation Coefficient (R)	0,9956

Serum: by applying the regression equation, the total systematic error estimated for decision levels of 80,0 µg/dL, 180,0 µg/dL and 270,0 µg/dL was -1.4%, 1.0% and 1.6%, respectively.

Precision:

Within-run precision was determined with 5 replicates for 10 days; between-run precision was determined in 3 condition, with 3 replicates for 5 days.

Samples (µg/dL)	Within-Run Precision		Samples (µg/dL)	Between-Run Precision	
	SD (µg/dL)	CV (%)		SD (µg/dL)	CV (%)
27.91	1.15	4.1%	28.93	1.42	4.9%
96.51	3.73	3.9%	94.60	3.09	3.3%
327.24	11.04	3.4%	324.56	12.88	4.0%

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Male	70 – 140 µg/dL
Female	80 – 155 µg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicatda.com.br.

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