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ENGLISH


INTENDED USE
Kit for the quantitative determination of the enzymatic activity of Alkaline Phosphatase in serum and plasma. Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Phosphorus Alkaline at alkaline pH hydrolyzes p-nitrophenylphosphate, releasing p-nitrophenol and phosphate. The catalytic concentration of the enzyme present in the sample is determined from the rate of formation of p-Nitrophenol, measured at 405nm.

SAMPLE – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and Plasma (heparin).
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.
Preservation: Alkaline Phosphatase in serum and plasma is stable for 7 days at a temperature of 4 to 8 °C and 2 months at a temperature of -20 °C.

PRODUCT DESCRIPTION

R 1	Diethanolamine ≥ 5.0 % v/v; Magnesium chloride ≥ 0.5 mmol / L; preservative.	<input checked="" type="checkbox"/>
R 2	Disodium PNPP ≥ 20 mmol / L; stabilizer; preservative.	<input checked="" type="checkbox"/>

Traceable to the specific absorbance of the chromogen p-nitrophenol using calibrated pipettes and handheld spectrophotometer that provide absolute values, based on modified methodology of the DGKC (German Society of Clinical Chemistry).

STABILITY IN USE

- The stability of the product (R1 and R2) in use is 18 months as long as followed by the recommended storage conditions (2 to 8 °C).
- The stability of working reagent is 28 days, as long as followed by the conditions of preparation and recommended storage (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Working Reagent (WR)
Mix in the ratio of: 4 parts of R1 + 1 part of R2. Homogenize gently.

B) OPERATING INTERVAL

The operating range of the product is 2.41 U/L to 700 U/L. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Control Internal Quality of Laboratory recommends the use of serum calibrator and control sera below:
Serum Calibrator - Autocal H 13.002.00
Normal Control Serum- Quantinorm **REF** 13.003.00
Serum Pathological Control - Quantialt 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Preheat the Working Reagent for three minutes at 37 °C.
2. Zero the equipment at 405nm with purified water.
3. Pipette into test tube:

	Volume
Work Reagent	1.0 mL
Sample	20 µL

4. Homogenize and insert into the thermostated cuvette at 37 °C. Set the stopwatch.
5. After 60 seconds, note the initial absorbance (A₀) and perform new readings every minute for 3 minutes (A₁, A₂ and A₃ respectively).

B) CALCULATIONS

Using the absorbance readings, calculate the mean absorbance change per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A_1 - A_0) + (A_2 - A_1) + (A_3 - A_2)}{3}$$

The catalytic activity of Alkaline Phosphatase in the sample is calculated by multiplying the ΔA / min by the corresponding factor:

$$\text{Alkaline Phosphatase (U/L)} = \Delta A/\text{min} \times 2720$$

Example:
A₀ = 1,179 A₁ = 1,225
A₂ = 1,266 A₃ = 1,312

$$\Delta A/\text{min} = \frac{(1,225 - 1,179) + (1,266 - 1,225) + (1,312 - 1,266)}{3}$$

ΔA/min = 0,044
Alkaline Phosphatase (U/L) = 0,044 x 2720 = 119,7 U/L

C) INTERPRETATION
Alkaline phosphatase activity is present in the body in membranes and cell surfaces located in the mucosa of the small intestine and in the proximal contorted tubules of the kidney, bone (osteoblasts), liver and placenta. Although its exact metabolic function is not understood, it appears that alkaline phosphatase is associated with the transport of lipids in the intestine and with the process of bone calcification. The enzyme is elevated

in different clinical conditions: cholestasis, diseases of the skeletal system (such as Paget's disease), hyperparathyroidism, rickets, osteomalacia, as well as cases of fractures and malignant tumors. Physiological increase may occur in children and young people as a result of increased osteoblastic activity in accelerated bone growth.

INTERFERING AND LIMITATIONS

Hemolysis, jaundice and Lipemia: Hemoglobin > 200 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 500 mg/dL.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 1.68 U / L / Limit of quantification: 2.41 U / L.

Analytical Specificity: The product specifically determines Alkaline Phosphatase in the presence of other interfering substances in the sample, up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 1,001x - 1,20$ and correlation coefficient $r = 0,9996$ was obtained. Using this equation the estimated total systematic error of - 0.50% for a level of 200 U/L and - 0.14% for a level of 500 U/L.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (U/L)	Repetitions	Intra-race Precision		Total Precision	
		SD (U/L)	%CV	SD (U/L)	%CV
174,103	80	0,865	0,5	0,916	0,5
371,608	80	0,626	0,2	1,329	0,4
633,232	80	0,704	0,1	0,704	0,1

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, calibrator, and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

	37 °C (U/L)
MEN	< 270 U/L
WOMEN	< 240 U/L

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): $\mu\text{kat} / \text{L}$
Alkaline Phosphatase (U / L) x 0.017 = Alkaline Phosphatase ($\mu\text{kat} / \text{L}$)

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer for reading at 405 nm.
- Water bath, thermostated at 37 ° C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Disposal, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels

of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of analysers.
The applications are available at www.biotechnicaltda.ind.br

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