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**ENGLISH**

**INTENDED USE**


Kit intended to determination of glucose in serum and plasma. Diagnostic use only.

**STORAGE AND HANDLING**

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

**WORKING PRINCIPLE**

The glucose present in the sample is oxidized by the enzyme glucose oxidase forming gluconic acid and hydrogen peroxide. It oxidizes the phenolic reactant in the presence of 4-aminoantipyrine, under the catalytic action of peroxidase, producing a pink compound (quinoneimine), which exhibits a maximum absorption at 505 nm. The color intensity is proportional to the glucose concentration in the sample.

**SAMPLE – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION**

**Sample Type:** Serum and plasma (fluoreto).

**Collection, handling and preparation:** Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Serum samples should be separated from the cell elements immediately.

**Preservation:** The glucose is stable in serum for 3 days at a temperature of 4 to 8 °C, 8 hours at a temperature of 20 to 25 °C and 24 hours at a temperature of -20 °C. Plasma glucose is stable for 3 days at a temperature of 20 to 25 °C.

**PRODUCT DESCRIPTION**

<b>R 1</b>	Phosphate buffer ≥ 40 mmol/L; 4-aminoantipyrine ≥ 0,2 mmol/L; phenol ≥ 1,0 mmol/L; Peroxidase ≥ 500 U/L; Glucose Oxidase ≥ 5000 U/L; activators; detergents and stabilizer.	<b>X</b>
<b>STD</b>	Glucose 100mg/dL; stabilizer and preservative. Traceable to NIST reference material 917c.	<b>X</b>

**STABILITY IN USE**

- The stability of the product (R1 e STD) in use is 18 months as long as followed by the recommended storage conditions (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

**TECHNICAL PROCEDURE**

**A) REAGENT PREPARATION**

**Reagent 1 (R1)**  
Ready-to-use Reagent

**B) OPERATING RANGE**

The product operating range is from 3,05 mg/dL to 400 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

**QUALITY CONTROL**

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:

Calibrator Serum - Autocal H	<b>REF</b>	13.002.00
Normal Control Serum – Quantinorm		13.003.00
Pathological Control Serum – QuantiAlt		13.004.00

**TEST PROCEDURE, CALCULATION AND INTERPRETATION**

**A) TEST PROCEDURE**

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 10 minutes at 37 °C.
3. Measure the absorbance of Standard and Sample against the blank at 505 nm. The final reaction is stable for 30 minutes.

**B) CALCULATIONS**

Glicose (mg/dL) =  $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{STD (mg/dL)}$

Example:  
STD concentration = 100 mg / dL  
Sample Absorbance = 0,219  
Absorbance of STD = 0,302  
Glicose (mg / dL) =  $\frac{0,219}{0,302} \times 100 = 72,5$  mg / dL

**With Calibration Factor (CF):**

CF =  $\frac{\text{STD (mg/dL)}}{\text{STD Absorbance}}$   
Glicose (mg/dL) = Sample Absorbance x CF  
Calibration Factor =  $\frac{100}{0,302} = 331$   
Glicose (mg/dL) = 0,219 x 331 = 72,5 mg/dL

**C) INTERPRETATION**

Glucose is the predominant carbohydrate in the peripheral blood, being the main source of cellular energy, stored under the forms of glycogen in the liver or fatty acids in adipose tissue. Their amount in the blood is controlled in narrow limits by several hormones, which have the function of reducing (insulin) or increasing (glucagon, adrenalina, cortisol and hormone of the growth). Hyperglycemia occurs in Diabetes mellitus caused by deficiency in the secretion or action of insulin. Other secondary factors can also elevate blood glucose levels, which include: pancreatitis, thyroid

dysfunction, renal failure or liver disease. Hypoglycemia is less prevalent, occurring in some diseases such as insulinoma, hypopituitarism or insulin-induced hypoglycemia.

#### INTERFERING AND LIMITATIONS

**Hemolyzed, Jaundice and Lipemic Sera:** Hemoglobin >500 mg/dL, Bilirubin > 15 mg/dL, Triglyceride > 1250 mg/dL interfere in the dosage.

**Drugs:** see recommended reference (Young, 2000).

#### PERFORMANCE CHARACTERISTICS

**Sensitivity:** Detection limit: 0,79 mg/dL / quantification limit: 3,05 mg/dL.

**Analytical Specificity:** The product determines glucose specifically in the presence of other interfering substances in the sample, up to the concentrations stated above

**Accuracy:** The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was  $y = 0,999x - 0,59$  and the correlation coefficient  $r = 0,9998$ . Using this equation, the total systematic error estimated is -0,84% to a level of 80U/L and -0,395% to a level of 200U/L.

**Precision:** Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained.

Samples (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
88,421	80	0,031	0,0	0,263	0,3
266,971	80	1,524	0,6	2,090	0,8
325,704	80	1,695	0,5	2,518	0,8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

#### RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

#### REFERENCE RANGES

Serum / Plasma (fasting)	Adult	74 - 99 mg/dL
	Child	60 - 100 mg/dL
	Premature infants	20 - 60 mg/dL
	Term Newborns	30 - 60 mg/dL

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): mmol/L Glucose (mg/dL) X 0,0556 = mmol/L

#### ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 505 nm (490 – 510 nm).
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

#### WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

#### QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br)

#### AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

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