HDL COLESTEROL

HDL CHOLESTEROL | HDL COLESTEROL

Ref. 10.005.00

FINALIDADE

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ENGLISH
INTENDED USE
Kit for the determination of HDL cholesterol in serum. Use in in vitro diagnostic.
STORAGE AND HANDLING
• Store at 2 to 8 °C.
• Protect from light.
• Stable until the kit expiration date that is printed on the label.
• Do not use reagents whose date has expired.

WORKING PRINCIPLE

The chylomicrons, the very low density lipoproteins (VLDL) and the low density lipoproteins (LDL) present in the sample are precipitated by phosphotungstate and magnesium ions. After centrifugation, cholesterol is determined in the supernatant containing HDL cholesterol (high density lipoproteins).

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.
Preservation: HDL cholesterol in serum is stable for 7 days at a temperature of 4-8 °C and stored at -20 °C for 2 months.

PRODUCT DESCRIPTION

R 1	Phosphotungstic acid > 1.0 mmol / L, Magnesium chloride > 10 mmol / L.	
STD	Butanol > 5% v / v, Triton X 100 5% v / v; stabilizer and concentration of HDL cholesterol indicated on the label.	X

STABILITY IN USE

- After opening, the product (R1 and STD) in use is stable up to the expiration date printed on the label, recommended storage conditions are followed (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1 and STD: Ready-to-use reagents.

B) OPERATING INTERVAL

The operating range of the product is from 3.0 mg / dL to 143 mg / dL. For higher values, dilute the supernatant with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration it is recommended to use the STD below:

HDL Cholesterol STD (included in kit)	REF	10.005.00
---------------------------------------	------------	-----------

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

Chylomicrons, VLDL and LDL precipitation:

1. Pipette in centrifuge tube:

Sample	250 µL
R1	250 µL

2. Shake vigorously for 30 seconds.
3. Centrifuge at 3,500 rpm for 15 minutes or until a clear supernatant is obtained.
4. Carefully collect the supernatant avoiding resuspension of the precipitate and carry out the dosage of cholesterol immediately.

HDL Cholesterol Dosage

Proceed to the determination of HDL cholesterol using Biotécnica R1 Cholesterol (CAT BT 10.004.00) as follows:

1. Pipette into test tubes

	Blank	STD	Sample
Purified water	100 µL	-	-
STD	-	100 µL	-
Supernatant	-	-	100 µL
Cholesterol R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize well and incubate the tubes for 10 minutes at 37 °C.
3. Measure the absorbance of the Standard and the Sample against the Blank at 505 nm. Color is stable for 20 minutes.

B) CALCULATIONS

To correct the dilution of the samples in the precipitation (1: 2), in the calculations, multiply the value of the Standard by 2.

HDL Cholesterol (mg/dL) = $\frac{\text{Sample Absorbance} \times [\text{Standard Concentration}(\text{mg/dL}) \times 2]}{\text{Standard Absorbance}}$

Exemplo:

Standard Concentration = 20 mg/dL
 Sample Absorbance = 0,244
 Standard Absorbance = 0,305
 HDL Cholesterol (mg/dL) = $\frac{0,244 \times [20 \times 2]}{0,305}$ = 32 mg/dL

With Calibration Factor:

With Calibration Factor = $\frac{[\text{Standard Concentration}(\text{mg/dL}) \times 2]}{\text{Standard Absorbance}}$
 HDL Cholesterol (mg/dL) = Sample Absorbance x Calibration Factor

Exemplo:

With Calibration Factor = $\frac{[20 \times 2]}{0,305}$ = 131
 HDL Cholesterol (mg/dL) = 0,244 x 131 = 32 mg/dL

C) INTERPRETATION

High-density lipoproteins (HDL) are responsible for the reverse transport of cholesterol from peripheral tissues to the liver. In the liver, cholesterol is converted to bile acids which are excreted in the intestine by bile. Monitoring of HDL cholesterol is of great clinical importance, as there is an inverse correlation between serum HDL concentration and the risk of atherosclerotic disease. Reduced concentrations of HDL cholesterol, along with high triglycerides, increase cardiovascular risk.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 150 mg / dL / Bilirubin > 15 mg / dL / Triglycerides > 1250 mg / dL, interfere with the dosage.
Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 0.828 mg / dL / Limit of quantification: 2,033 mg / dL.
Analytical Specificity: The product specifically determines HDL cholesterol in the presence of other interfering substances in the sample up to the concentrations reported above.
Accuracy: The method was compared with a similar method by the determination of 30 samples in duplicate. The regression equation Y = 1.031x - 0.490 and correlation coefficient r = 0.9901 was obtained. Using this equation the estimated total systematic error of 1.47% for a level of 30.0 mg / dL and 2.61% for a level of 100.0 mg / dL.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
10,6	80	0,107	1,0	0,223	2,1
50,9	80	0,316	0,6	1,317	2,6
127	80	0,515	0,4	0,695	0,5

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Age Group	Desirable reference values (mg / dL)
Grown up (> 20 anos)	> 40
Children and adolescents	> 45

NOTE: Reference values are desirable according to the Brazilian Consensus for the Normalization of Laboratory Determination of the Lipid Profile.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 505 nm (490 -540).
- Water bath, thermostated at 37 ° C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaitda.com.br

AUTOMATION

The HDL Cholesterol dosing procedure is automated in most analyzers. The applications are available at www.biotecnicaitda.ind.br

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- CASTELLI WP et al. HDL Cholesterol and other lipids in coronary heart disease. *Circulation ahajournals.org*. 55;767-772 (1977).
- GORDON T et al., High density lipoprotein as a protective factor against coronary heart disease. *Am. J. Med.*, 62: 707-714 (1977).
- FALUDI AA, IZAR MCO, SARAIVA JFK, CHACRA APM, BIANCO HT, AFIUNE NETO A et al. Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose -2017. *Arq Bras Cardiol* 2017; 109(2Supl.1):1-76.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D. S. **Effects of drugs on clinical laboratory tests - vol. 2, 5 ed.** Washington DC: AAC Press, 2000.
- WARNICK, GR; WOOD, PD. National cholesterol education program recommendations for measurement of high-density lipoprotein cholesterol: executive summary. *Clin. Chem.* V.41. 10, 1427-1433 (1995).
- NCEP- Detection, evaluation and treatment of high blood cholesterol in adults (Adult Treatment Panel III). NHH Publication 02-5215, Bethesda, MD, 2002.
- Consenso Brasileiro para a Normalização da Determinação Laboratorial do Perfil Lipídico, 2016.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABELA DE SÍMBOLOS INTERNACIONALES

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