IBC
IBC | UIBC
Ref. 12.010.00

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ENGLISH
INTENDED USE
Kit for the determination of iron binding capacity in serum. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
The iron ions present in R1 saturate the available transferrin sites in the sample. The excess of unbound iron reacts with the chromogen ferrozine to form a colored complex, which exhibits a maximum absorption at 560 nm. The difference between the total amount of iron contained in R1 and the quantity determined for the sample corresponds to the Latent Iron Binding Capacity of the sample. The Total Iron Binding Capacity is the sum of the Latent Iron Binding Capacity and the Serum Iron.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Serum.
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Centrifuge and separate the serum within two hours after collection.
Preservation: The sample is stable for 6 months at - 20 °C.

PRODUCT DESCRIPTION

R 1	TRIS Buffer ≥ 100 mmol / L; ferric chloride ≥ 10 µmol / L; detergent; stabilizer; preservative.	
R 2	Tampão TRIS ≥ 10 mmol/L; Ácido Ascórbico ≥ 20 mmol/L; Ferrozine ≥ 2,0 mmol/L; estabilizante.	
CAL	Lyothesized human serum containing human transferrin, stabilizers and preservatives. Standard value printed on the bottle label. Traceable to NIST reference material SRM-937.	

STABILITY IN USE

- **R1 and R2:** After opening, the product in use is stable to the expiration date printed on the label, recommended storage conditions are followed (2 to 8 °C).
- **CAL:** After reconstitution, the calibrator is stable for 5 days if stored at 2 to 8 °C and 30 days at - 20 °C.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT
A) REAGENT PREPARATION
R1 and R2: Ready-to-use reagents.
CAL: Freeze-dried. Reconstitute with 1.5 mL of purified water. Leave to stand for 30 minutes. Mix by gentle inversion avoiding foaming.

B) OPERATING INTERVAL
The operating range of the product is 39.8 µg / dL at 469 µg / dL.
For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL
The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For the Internal Quality Control Laboratory it is recommended to use the control sera below:

Normal Control Serum- Quantinorm	REF	13.003.00
Serum Pathological Control - Quantialt		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE
1 - Pipette into test tube.

	Blank	Test	Calibrator
R1	0,75 mL	0,75 mL	0,75 mL
Sample	----	0,1 mL	----
CAL	----	----	0,1 mL
Purified Water	0,1 mL	----	----

2. Homogenize gently. Incubate for 5 minutes at 37 °C.
3. Measure the absorbance of the Blank, Test and Calibrator at 560nm (525-575), setting the zero with purified water. This will be the A1 absorbance.
4. Add:

	Blank	Test	Calibrator
R2	0,2 mL	0,2 mL	0,2 mL

5. Homogenize gently and incubate 5 minutes at 37 °C.
6. Measure the Absorbance of Blank, Test and Calibrator at 560nm (525-575). This will be the absorbance A2. The color is stable for 30 minutes.

B) CALCULATIONS
The A1 absorbances must be corrected to the final volume of the reaction, obtaining A1COR, according to the following equation:

$$A_{1COR} = A_1 \times 0,809$$

Determine the ΔA of the Blank, Test and CAL: A₂ - A_{1COR}.
For calculations use the absorbance difference (A₃) between CAL / Test and White:

A₃ CAL: ΔA White - ΔA CAL
A₃ Test: ΔA White - ΔA Test

Latent Iron Binding Capacity:

$$CLLF (\mu\text{g} / \text{dL}) = \frac{A_3 \text{ Test}}{A_3 \text{ CAL}} \times \text{CAL Concentration}$$

Example:
White:
A₁ = 0.003 A_{1COR} = 0.003 x 0.809 = 0.002
A₂ = 0.700 ΔA White = 0.700 - 0.002 = 0.698

CAL:
CAL concentration = 287 µg/dL
A₁ = 0,056 A_{1COR} = 0,056 x 0,809 = 0,045
A₂ = 0,645 ΔA CAL = 0,645 - 0,045 = 0,600
A₃ = 0,698 - 0,600 = 0,098

Test:
A₁ = 0,027 A_{1COR} = 0,027 x 0,809 = 0,022
A₂ = 0,658 ΔA Teste = 0,658 - 0,022 = 0,636
A₃ = 0,698 - 0,636 = 0,062

CLLF (µg/dL) = $\frac{0,062}{0,098} \times 287 = 182 \mu\text{g/dL}$

Total Iron Bonding Capacity:
CTLF ($\mu\text{g} / \text{dL}$) = Serum Iron + CLLF

Transferrin Saturation%:

$$\text{Transferrin Saturation (\%)} = \frac{\text{Serum Iron}}{\text{CTLF}} \times 100$$

C) INTERPRETATION

About 1 mg of iron is absorbed daily. Once in the mucosal cells, Fe²⁺ ions bind to transporters. Before passing into the plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bind to apotransferrin. The apotransferrin-Fe³⁺ complex is called transferrin which carries the iron ions in the blood plasma. Approximately one third of the transferrin binding sites are occupied by Fe³⁺, so the additional amount of iron that can be bound is latent iron binding capacity (CLLF). The sum of serum iron with CLLF represents the total iron binding capacity (CTLF). CTLF varies in iron metabolism disorders. In iron deficiency anemia CTLF is elevated and transferrin saturation is low. Low serum iron concentrations associated with low CTLF may occur in cases of chronic anemia, malignant tumors or infections.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 150 mg / dL / Bilirubin > 20 mg / dL / Triglycerides > 1000 mg / dL, interfere with the assay.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 17.61 $\mu\text{g} / \text{dL}$ / Limit of quantification: 39.77 $\mu\text{g} / \text{dL}$.
Analytical Specificity: The product specifically determines iron binding capacity in the presence of other interfering substances in the sample up to the concentrations reported above.
Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 1.004x - 1.04$ and correlation coefficient $r = 0.9945$ was obtained. Using this equation the estimated total systematic error of -0.293% for a level of 150.00 $\mu\text{g} / \text{dL}$ and 0.053%, to a level of 300.00 $\mu\text{g} / \text{dL}$.
Accuracy: it was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples ($\mu\text{g}/\text{dL}$)	Repetitions	In-run Precision		Total Precision	
		SD ($\mu\text{g}/\text{dL}$)	%CV	SD ($\mu\text{g}/\text{dL}$)	%CV
105,5	80	1,16	1,1	3,25	3,1
173,8	80	1,16	0,7	4,45	2,6
336,2	80	2,14	0,6	4,32	1,3

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Latent Iron Capacity - CLLF	110 to 370 $\mu\text{g}/\text{dL}$
Total Iron Capacity - CTF	250 to 450 $\mu\text{g}/\text{dL}$
Transferrin Saturation index - IST	15 to 55 %

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): $\mu\text{mol} / \text{L}$

CTLF ($\mu\text{g} / \text{dL}$) $\times 0,179 =$ CTLF ($\mu\text{mol} / \text{L}$)

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 560 nm (525-575).
- Water bath, thermostated at 37 ° C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use all BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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