LDL Colesterol

LDL Cholesterol | LDL Colesterol

Ref. 10.015.00

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ENGLISH

INTENDED USE

Kit for the determination of LDL cholesterol in serum and plasma. Use in *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

R1 contains polyvinyl sulfonic acid (PVS) and polyethylene glycol methyl ether (PEGME) which forms a complex with LDL-VLDL-chylomicrons, leaving HDL cholesterol free. This is consumed by Cholesterol oxidase and Cholesterol esterase. After the addition of R2 the LDL is released from the PVS / PEGME complex by specific detergents. Free LDL reacts with the enzymes producing hydrogen peroxide which is quantified by the Trinder reaction. The intensity of the color produced is proportional to the amount of LDL cholesterol present in the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum or plasma (heparin and EDTA).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: LDL in serum and plasma is stable for 7 days at a temperature of 4-8 °C and 2 months at -20 °C.

PRODUCT DESCRIPTION

R 1	MES buffer pH 6.5; polyvinylsulfonic acid, polyethylene glycol methyl ester, Magnesium Chloride, Detergents, EDTA, 4-aminoantipyrine, Cholesterol esterase, Cholesterol Oxidase and Peroxidase.
	MES buffer pH 6.5, EDTA, TODB N Detergent, N-bis (4-sulfobutyl) -3-metilaniilina
R 2	

Human serum, preservatives. Value printed on the bottle label. Traceable to CDC reference method for LDL Cholesterol.

STABILITY IN USE

- R1 and R2: Once opened, reagents in use are stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 ° C) are followed.
- CAL: After reconstitution, CAL is stable for 14 days if stored at 2 to 8 ° C and 3 months at -20 ° C.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1 and R2: Ready-to-use reagents.

CAL: Freeze-dried. Reconstitute with 1 mL of purified water. Leave to stand for 30 minutes. Homogenize gently by inversion avoiding foaming.

B) OPERATING INTERVAL

The operating range of the product is from 5,5 mg / dL to 241 mg / dL.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Laboratory Quality Control it is recommended to use the calibrator serum and control sera below:

LDL Cholesterol CAL (included in the kit)	10.015.00
Normal Control Serum - Quantinorm	13.003.00
Serum Pathological Control - Quantialt	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

REACTION PARAMETERS	
Reaction Types	Final/Growing Point
Wave length	600 nm
Temperature	37 °C
Purified Water (white) / sample/CAL*	3 µL
R1*	225 µL
Reading 1 (A ₁)	5 minutes
R2*	75 µL
Reading 2 (A ₂)	5 minutes
Calibration	Linear

* The volumes may be modified provided that the established ratio is maintained.

B) CALCULATIONS

LDL Cholesterol (mg / dL) = $\frac{\Delta A \text{ Sample} - \Delta A \text{ White}}{\Delta A \text{ Calibrator} - \Delta A \text{ White}}$ x CAL Value (mg / dL)

Where: $\Delta A = A_2 - A_1$

C) INTERPRETATION

Very low density lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on VLDL (Very Low Density Lipoproteins). The elimination of LDL cholesterol from plasma occurs primarily by the liver in parenchymal cells through specific LDL receptors. Most of the cholesterol present in atherosclerotic plaques originates from LDL. For this reason LDL cholesterol is considered an important clinical predictor of coronary atherosclerosis.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 1000 mg / dL / Bilirubin > 40 mg / dL / Triglycerides > 1500 mg / dL, interfere in the dosage.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Limit of detection: 1,501 mg / dL / Limit of quantification: 5,474 mg / dL.

Analytical Specificity: the product specifically determines LDL cholesterol in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 37 samples in duplicate. The regression equation $Y = 0,994x + 2,02$ and correlation coefficient $r = 0,9953$ was obtained. Using this equation the total systematic error estimated was 3.28% for a level of 52 mg / dL and 0.32% for a level of 220 mg / dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	In-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
32,401	80	0,580	1,8	1,019	3,1
62,944	80	0,689	1,1	1,023	1,6
160,263	80	0,947	0,6	1,423	0,9

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES (THERAPEUTICAL TARGET VALUES)

Suggested therapeutic target values for cardiovascular risk category estimated by the Doctor for adults over 20 years:

RISK	THERAPEUTIC TARGET (mg/dL)
Low	< 130
Intermediary	< 100
High	< 70
Very High	< 50

NOTE: According to the Brazilian Consensus for the Normalization of Laboratory Determination of the Lipid Profile, reference values for LDL cholesterol were replaced by the values of therapeutic target, stratified according to the calculated cardiovascular risk.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 600 nm .
- Glass and / or automatic pipettes.
- Clock or Stopwatch.

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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- BADIMON JJ, BADIMON L, FUESTER V. Regression of atherosclerotic lesions by high density lipoprotein plasma fraction in the cholesterol-fed rabbit. *Journal of Clinical Investigation*, 85: 1234-41 (1990).
- BARR DP, RUSS EM, EDER HA. Protein-lipid relationships in human plasma. *Am. J. Med.*, 11:480-493 (1951).
- CASTELLI WP et al. HDL Cholesterol and other lipids in coronary heart disease. *Circulation ahajournals.org*. 55;767-772 (1977).
- CROUSE JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease. *J. Lipid Res.*, 26, 566-574 (1985).
- FALUDI AA, IZAR MCO, SARAIVA JK, CHACRA APM, BIANCO HT, AFIUNE NETO A et al. Atualização da Diretriz Brasileira de Dislipidemias e Prevenção da Aterosclerose - 2017. *Arq Bras Cardiol* 2017; 109(2Supl.1):1-76.
- GORDON T et al., High density lipoprotein as a protective factor against coronary heart disease. *Am. J. Med.*, 62, 707-714 (1977).
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACC Press, 2000.

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