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ENGLISH
INTENDED USE
Kit intended to determination of lactate in plasma. In vitro Diagnostic use only.
STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
The lactate in the sample is oxidized by the enzyme lactate oxidase, and in the presence of oxygen, produces pyruvate and hydrogen peroxide. This, in the presence of a phenolic reagent (TOOS) and 4-aminoantipyrine and peroxidase produces a chromogen with maximum absorption at 546 nm.
Lactate + O₂ $\xrightarrow{\text{Lactate Oxidase}}$ Pyruvate + H₂O

H2O2 + 4-aminoantipyrine + TOOS $\xrightarrow{\text{Peroxidase}}$ Chromogen + 2 H2O

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Plasma (fluoride).
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practices. The samples to be analyzed should be treated as potentially infectious material.
Preservation: The sample should be separated in up to 15 minutes. The plasma is stable for 8 hours at 15-25 ° C, 14 days if stored at 2-8 ° C or 30 days at -20 °C.

PRODUCT DESCRIPTION

R 1 Buffer Pipes ≥ 10 mmol / L; Lactate oxidase ≥ 100 U/L; Peroxidase ≥ 500 U/L; 4-aminoantipyrine ≥ 0.1 mmol/L; TOOS ≥ 0.5 mmol / L; Cofactor; Stabilizers; Activators; detergent; Preservative.
Phosphate Buffer ≥ 20 mmol/L; Bovine albumin ≥ 5% w/v; Stabilizer; Preservative and Lactate at a concentration equivalent to 40 mg/dL. Traceable to NIST SRM 84L reference material and IC Assay NIST SRM

CAL Phosphate Buffer ≥ 20 mmol / L; Bovine albumin ≥ 5% w/ v; Stabilizer; Preservative and lactate at the specified concentration on the vial label. Traceable to NIST SRM 84L reference material and IC Assay NIST SRM.

CONTROL Phosphate Buffer ≥ 20 mmol / L; Bovine albumin ≥ 5% w/ v; Stabilizer; Preservative and lactate at the specified concentration on the vial label. Traceable to NIST SRM 84L reference material and IC Assay NIST SRM.

STABILITY IN USE

- The product stability (R1, CAL and CONTROL) in use is 14 months, since recommended storage conditions are followed (2 to 8 ° C)
- Reagents should remain outside the specified temperature only for the time required to perform the tests

PRODUCT TREATMENT AND HANDLING
A) REAGENT PREPARATION
Reagent 1, Calibrator and Control (R1 / CAL / CONTROL)
Ready for use.
B) OPERATING INTERVAL
The product operating range is from 2.32 mg/dL to 120.00 mg/dL. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL
The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and the control sera below:
Serum Calibrator - Autocal H 13.002.00
Normal Control Serum - Quantinorm 13.003.00
Serum Pathological Control - Quantialt 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE
1. Pipette in the assay tubes:

	Blank	CAL	Sample
CAL	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize well and incubate the tubes for 5 minutes at 37 ° C.
3. Measure the Absorbance of Calibrator (Ac) and Sample (Aa) versus White at 546 nm.
B) CALCULATIONS
Calculation for Plasma:
Lactate (mg / dL) = Sample Absorbance x Concentration Calibrator (mg/dL) / Calibrator Absorbance

Example:
Calibrator Concentration = 40 mg/dL
Sample Absorbance = 0.225
Calibrator Absorbance = 0.603
Lactate (mg/dL) = $\frac{0,225}{0,603} \times 40 = 14,9$ mg/dL
With Calibration Factor:
Calibration Factor = $\frac{\text{Calibrator Concentration (mg/dL)}}{\text{Calibrator Absorbance}}$
Lactate (mg / dL) = Sample Absorbance x Calibration Factor

Example:
Calibration Factor = $\frac{40}{0,603} = 66,3$
Lactate (mg/dL) = $0,225 \times 66,3 = 14,9$ mg/dL.
C) INTERPRETATION
Lactate is an intermediate product in the metabolism of carbohydrates that predominantly accumulates in the blood when tissue oxygenation decreases. It

originates mainly from the skeletal muscle, brain, skin, kidney marrow and erythrocytes. The concentration of blood lactate depends on the rate of production in these tissues and the metabolic rate in the liver and kidneys. The liver uses approximately 65% (75 g / day) of total basal lactate produced predominantly in gluconeogenesis. Anaerobic glycolysis markedly increases lactate in the blood, especially in prolonged exercise. The anoxia resulting from shock, pneumonia and congestive heart failure; Thiamine deficiency and diabetic ketoacidosis, among other causes, cause increased lactate and blood pyruvate. High blood lactate concentration results in lactic acidosis, which can occur in renal failure and leukemia.

INTERFERES OR LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 500 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 1000 mg/dL.

Medications: consult recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 0.84 mg/dL / Limit of quantification: 2.32 mg/dL.

Analytical Specificity: The product specifically determines lactate in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0.990x - 0.16$ and correlation coefficient $r = 0.9983$ was obtained. Using this equation the estimated total systematic error of -1.32% for a level of 50 mg/dL and -1.16% for a level of 100 mg/dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mg/dL)	Repetition s	Within-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
25,281	80	0,244	1,0	0,568	2,2
35,246	80	0,247	0,7	0,524	1,5
105,411	80	0,230	0,2	0,631	0,6

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The water level of the water bath should be higher than the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Sample	mg/dL	mmol/L	Specification
Plasma (Fluoride)	4,5-19,8	0,5-2,2	Venous
	4,5-14,4	0,5-1,6	Arterial

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI):

Lactate (mg / dL) x 0.111 = Lactate (mmol/L)

MATERIAL NEEDED TO PERFORM THE TEST

- Spectrophotometer or photometer for reading at 546 nm (530 - 570).
- Water bath, thermostated at 37 ° C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Assay tubes.

WARNINGS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of the leftovers of the reactions according to the Good Clinical Laboratory Practices (GCLP) and Health Service Waste Management Program (HSWMP).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotécnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@Biotechnicaltda.com.br

AUTOMATION

This procedure is compatible to the most types of analysers. The applications are available at www.biotechnicaltda.ind.br

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- HUTCHESON, A. et al. Measurement of lactate in cerebrospinal fluid in investigation of inherited metabolic disease. **Clin. Chem.** v.43, n.1, p.158-161, 1997.
- GUYTON, A. C.; HALL, J. E. **Fisiologia humana e mecanismos das doenças**. 6ª ed. Rio de Janeiro: Guanabara Koogan. 1998. 580p.
- FOSSATI, P.; PRENCIPE, L.; BERT G. Use of 3,5-dichloro-2-hydroxibenzensulfonic acid/4-aminophenazone chromogenic system in direct enzymic assay of uric acid in serum and urine. **Clin. Chem.** v.26, p.227-231, 1980.
- BARHAM D.; TRINDER P. An improved colour reagent for the determination of blood glucose by oxidase system. **Analyst** v.27, p.142-145, 1972.
- SHINOJO, N. et al. Test strip method for measuring lactate in whole blood. **Clin. Chem.** v.35, n.9, p.1992-1994, 1989.
- BURTIS, C. A.; ASHWOOD, E. R.; BRUNS, D. E. **Tietz Fundamentos de Química Clínica**, Saunders Elsevier, 6 ed, 2008.
- BURTIS, C. A.; ASHWOOD, E. R. **Tietz Fundamentos de Química Clínica**, Guanabara Koogan, 4 ed, 1998.
- YOUNG, D. S. **Effects of drugs on clinical laboratory tests** - vol. 2, 5 ed. Washington DC: AACCC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. **Clin. Chem.** v.27 p.493-501, 1981.
- WESTGARD, James O.; LAHMEYER, Brenda L.; BIRNBAUM, Marvin L. Use of the Du Pont "Automatic Clinical Analyzer" in Direct Determination of Lactic Acid in Plasma Stabilized with Sodium Fluoride. **Clin. Chem.**, Cincinnati, v. 18, nº. 11, p.1334-1338, nov. 1972.

TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS / TABLA DE SÍMBOLOS INTERNACIONALES			
	Consultar Instruções de Uso Consult instructions for use Consultar Instrucciones de Uso	Σ	Conteúdo suficiente para <-> testes Contains sufficient for <->tests Contenido suficiente para <->ensayos
REF	Código Code Código		Conteúdo suficiente para <-> testes Contains sufficient for <->tests Contenido suficiente para <->ensayos
LOT	Número de lote Batchcode Denominación de lote		Límite de temperatura Temperature limitation Temperatura límite
IVD	Para uso diagnóstico <i>in vitro</i> For <i>in vitro</i> diagnostic medical device Para uso em diagnóstico <i>in vitro</i>		Novo / Irritante Harmful / Irritant Novo / Irritante
R <->	Reagente e seu número/abreviação Reagent and its number/abbreviation Reactivo y su número/abreviación		Risco biológico Biological risk Riesgo biológico