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SENSIBLE PROTEIN | PROTEÍNA SENSÍVEL
Ref. 10.016.00 (KIT) / 10.016.00 (KIT COM CONTROLE)

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ENGLISH

INTENDED USE

Kit intended to determine protein in urine and liquor samples. Use in vitro diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagent ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

WORKING PRINCIPLE

Method: Pyrogallol Red

The protein reacts with pyrogallol red and sodium molybdate, in acidic medium, forming a colored complex that can be spectrophotometrically measured at 600 nm. The intensity of the color is proportional to the concentration of protein in the sample.

SAMPLE: TYPE, COLLECTION, HANDLING, PREPARATION AND STABILITY

Sample Type: urine and cerebrospinal fluid.

Collection and Handling: collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

Preservation:

	Temperature	Stability Period
Urine	4 to 8 °C	7 days
	-20 °C	1 month
Cerebrospinal fluid	4 to 8 °C	6 days
	-20 °C	> 1 year

Preparation:

Urina: homogeneize the sample and separate 20 mL aliquot. Centrifuge for 10 minutes at 3000 rpm and utilize the supernatant to proceed with the assay.
Cerebrospinal Fluid: CSF samples should be centrifuged.

PRODUCT DESCRIPTION

R 1	Succinate buffer > 10 mmol / L; sodium molybdate > 10 µmol / L; pyrogallol red > 10 µmol / L; preservative; detergent; stabilizers. Pipes buffer; preservatives; stabilizer; Bovine albumin at a concentration corresponding to 50,0 mg / dL protein. Traceable to NIST reference material 927d.	
STD	Phosphate buffer; stabilizers; detergent; preservatives and bovine albumin. (INCLUDED IN THE KIT WITH CONTROL)	
CONTROL		

QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the control below:

SENSIBLE PROTEIN CONTROL
(INCLUDED IN THE KIT WITH CONTROL) **REF** 10.016.00

NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 600 nm (580 – 620 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Sample	STD
STD	---	---	50 µL
Sample	---	50 µL	---
R1	1,0 mL	1,0 mL	1,0 mL

- Homogenize and incubate for 10 minutes at 37 °C. Remove it from the water bath and let it rest in temperature of 15 to 30 °C before reading.
- Measure the Standard's and Sample's absorbance against the blank at 600 nm (580 – 620 nm). The final reaction is stable for 15 minutes.

B) CALCULATIONS

Protein (mg/dL) = $\frac{\text{Sample Absorbance} - \text{Blank}}{\text{STD Absorbance} - \text{Blank}} \times \text{STD Concentration (mg/dL)}$

Calculations with the Calibration Factor:

Calibration Factor = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance} - \text{Blank}}$

Protein (mg/dL) = (Sample Absorbance - Blank) x Calibration Factor

Calculations for urine samples:

Protein (mg/24 hours) = $\frac{\text{Protein (mg/dL)} \times \text{Urinary Volume (mL)}}{100}$

Automation: this product is compatible to most types of automatic analyzers. Instrument settings are available at www.biotecnicaldta.ind.br

C) INTERPRETATION

Proteinuria is defined as the presence of excess protein in the urine, which occurs in many types of kidney disease, after intense exercise or dehydration. The causes may be due to increased glomerular permeability, altered tubular reabsorption, elevated low-molecular-weight plasma proteins, abnormal secretion of proteins in the urinary tract and post-renal bleeding.

In cerebrospinal fluid, the quantification of total proteins is used mainly for the detection of diseases in the Central Nervous System (CNS) related to increased permeability of the blood-brain barrier. High values can be found in infections, intracranial hemorrhage, multiple sclerosis, Guillain-Barré syndrome, malignancies, some endocrine abnormalities, use of some medications and in various inflammations. Extremely low values are found when there is fluid loss in the CNS.

PERFORMANCE CHARACTERISTICS

Operating range
4,80 a 97,00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0,9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0,50 mg/dL	4,80 mg/dL

Analytical Specificity	
Hemoglobin	Bilirubin
Any concentration	> 1,5 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 0.992x + 0.028
Correlation Coefficient (R)	0.9975

Urina: by applying the regression equation, the total systematic error estimated for decision levels of 30.0 mg/dL and 70.0 mg/dL was -0.71% and -0.76%, respectively.

Precision:

Determined with two runs in duplicate per day for 20 days.

Samples (mg/dL)	Replicates	Within-Run Precision		Total Precision	
		SD (mg/dL)	CV (%)	SD (mg/dL)	CV (%)
15.20	80	0.276	1.80	0.752	5.00
47.50	80	0.518	1.10	1.905	4.00
87.10	80	0.693	0.80	1.730	2.00

CV: Coefficient of variation; SD: Standard deviation.

RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

REFERENCE RANGES

Urine		Cerebrospinal Fluid	
Adults	1 - 14 mg/dL	0 to 30 days	15 - 153 mg/dL
	<100 mg/day	3 to 6 months	15 - 48 mg/dL
Pregnant	< 150 mg/day	Up to 2 years	15 - 50 mg/dL
		Older than 2 years	15 - 45 mg/dL

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

Conversion to International System of Units (g/L):

Protein (mg/dL) x 0,01 = Protein (g/L)

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to sac@biotecnicatda.com.br.

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- BIRMINGHAM, D. J. et al. Relationship between Albuminuria and Total Proteinuria in Systemic Lupus Erythematosus Nephritis: Diagnostic and Therapeutic Implications. *Clinical Journal Of The American Society Of Nephrology*, [s.l.], v. 3, n. 4, p.1028-1033, 9 abr. 2008
- BURTS, C. A.; ASHWOOD, E. R. **Tietz textbook of clinical chemistry**, 4 ed. Philadelphia: W.B.Saunders, 1994. 836p.
- GNUTZMANN, L. V.; PLEWKA, J.; SULDOFSKI, M. T.; FELISBERTO, M.; NESI, V. Análise dos valores de referência do líquido cefalorraquídeo. *RBAC*, v. 48, n. 3, p. 189-197, 2016.
- MARSHALL, Thomas; WILLIAMS, Katherine M.. Total Protein Determination in Urine: Elimination of a Differential Response between the Coomassie Blue and Pyrogallol Red Protein Dye-binding Assays. *Clinical Chemistry*, v. 46, n. 3, p.392-398, 2000.
- SOLDIN, J. O.; WONG, E. C.; BRUGNARA, C. Pediatric Reference Intervals, 7ª.edition Washington: AACCPress, 2011.p.193.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

- YALAMATI, Padma; KARRA, Madhu Latha; BHONGIR, Aparna V.. Comparison of Urinary Total Proteins by Four Different Methods. *Indian Journal Of Clinical Biochemistry*, [s.l.], v. 31, n. 4, p.463-467, 8 fev. 2016. Springer Nature.
- YOUNG, D.S. **Effects of drugs on clinical laboratory tests - vol. 2**, 5 ed. Washington DC: AACCPress, 2000.

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