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ENGLISH

INTENDED USE



Kit intended to determination of AST/GOT activity in serum and plasma (EDTA and heparina). Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The Aspartate Transaminase (AST or GOT) catalyses the transference of amine group from aspartate to α -ketoglutarate. This transference produces oxaloacetate and glutamate. The activity is then determined by coupling this reaction to Malate Dehydrogenase (MDH) enzyme that transforms the oxaloacetate product in malate consuming NADH, which is measured at 340 nm.

SAMPLE – PREPARATION AND STABILITY

Sample Type: Serum and plasma (EDTA and heparina)

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The AST-GOT in serum and plasma is stable for 7 days if stored at a temperature of 4 to 8 °C.

PRODUCT DESCRIPTION

Tris buffer ≥ 80 mM, L-aspartic acid ≥ 100 mM, alpha-ketoglutarate ≥ 10 mM, Malate Dehydrogenase ≥ 100 U/L, Lactate Dehydrogenase ≥ 500 U/L, activators, stabilizers, detergentes, preservative.

NADH ≥ 0,5 mM, stabilizer, preservative.

The determination of AST-GOT is traceable to reference material ERM – AD457/IFCC.

STABILITY IN USE

- The stability of the product (R1 and R2) in use is 18 months as long as followed by the recommended storage conditions (2 to 8 °C).
- The stability of working reagent is 14 days, as long as followed by the conditions of preparation and recommended storage (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Work Reagent (WR)

Mix in proportion: 4 parts of R1 + 1 part of R2.

Homogenize it gently. The Work Reagent is stable 14 days at 2 to 8 °C.

B) OPERATING RANGE

The product operating range is from 3,05 U / L to 440 U / L. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:

Calibrator Serum - Autocal H

Normal Control Serum – Quantinorm

Pathological Control Serum – QuantiAlt

REF	13.002.00
REF	13.003.00
REF	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Heat the reagent for 3 minutes at 37 °C.

2. Pipette in the assay tube:

	Volume
Work Reagent	1.0 mL
Sample	100 µL

3. Homogenize and insert it immediately in the thermostated cuvette at 37 °C.
4. After 1 min. note the initial absorbance (A₀) and read again after exactly 1, 2 and 3 minutes (A₁, A₂ and A₃, respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min} = \frac{(A_0 - A_1) + (A_1 - A_2) + (A_2 - A_3)}{3}$$

The activity of the GOT is then calculate by multiplying the ΔA/min by the factor below:

$$\text{GOT (U/L)} = \Delta A/\text{min} \times 1746$$

Example:

A₀ = 1,268; A₁ = 1,228; A₂ = 1,189; A₃ = 1,152

$$(\Delta A/\text{min}) = \frac{(1,268 - 1,228) + (1,228 - 1,189) + (1,189 - 1,152)}{3}$$

GOT (U/L) = 0.039 x 1746 = 68,1 U/L

C) INTERPRETATION

The Glutamate Oxaloacetate Transaminase (GOT) or Aspartate Transaminase (AST) is an intracellular enzyme present in high quantities in cardiac and skeletal muscles, liver, kidney and brain. Elevated levels of GOT help the diagnostic of cardiac, hepatic and muscle diseases. The activity of this enzyme in the myocardial infarction is elevated within the first 12 hours, reaching a peak in 24 hours and come back to normal in the fifth day. Discrete elevations are observed in pregnancy. Elevated levels are also found in hepatic necrosis, hemolytic anemias, acute pancreatitis, cirrhosis, hepatitis, obstructive jaundice, mononucleosis, hypothyroidism, brain necrosis and trauma, severe burns, muscle dystrophies, skeletal muscle injuries, cardiac caterism and angioplasty. Several common drugs may increase the level of GOT (isoniazide, erythromycin, progesterone, anabolic steroids, etc.), being then used in monitoring therapies which uses hepatotoxic drugs.

INTERFERING AND LIMITATIONS

Anticoagulants: citrate, sodium fluoride and sodium oxalate may interfere.

Hemolyzed, Jaundice and Lipemic Sera: Hemoglobin > 200 mg/dL, Bilirubin > 40 mg/dL, Triglyceride > 350 mg/dL.

Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 2,41 U / L / quantification limit: 3,05 U / L.

Analytical Specificity: The product determines AST-GOT specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 1,006x + 0,53$ and the correlation coefficient $r = 0,9996$. Using this equation the total systematic error estimated is 1,66% to a level of 50 U/L and 0,95% to a level of 150 U/L.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained.

Samples (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
64,17	80	1,05	1,6	1,18	1,8
240,49	80	1,27	0,5	1,27	0,5
424,86	80	2,77	0,7	3,22	0,8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

	ACTIVITY AT 37 °C
MEN	< 35 U/L
WOMEN	< 31 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): GOT (U/L) X 0,017 = GOT (μKat/L)

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 340 nm.
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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- SCHUMANN, G.; et al. IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37 °C. Part 5. Reference Procedure for the Measurement of Catalytic Concentration of Aspartate Amino-transferase. *Clin. Chem. Lab. Med.* v. 40 p.725-733, 2002.
- BURTIS, C. A.; ASHWOOD, E. R. Tietz textbook of clinical chemistry, 4 ed. Philadelphia: W.B.Saunders, 1994. 836p.
- YOUNG, D. S. **Effects of drugs on clinical laboratory tests - vol. 2**, 5 ed. Washington DC: AACC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

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