ALBUMINA
ALBUMIN / ALBÚMINA
Ref. 10.002.00

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ENGLISH

INTENDED USE

Kit intended to determination of albumin in serum. Diagnostic use only.

STORAGE AND HANDLING

- Store at 15 to 30 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Albumin reacts with bromocresol green at acid pH to form a blue-green complex that is measured at 630 nm. The intensity of the color is proportional to the concentration of albumin in the sample.

SAMPLE – PREPARATION AND STABILITY

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Serum albumin is stable for 3 days if stored at 4 to 8 ° C and 7 days if stored at -20 ° C

PRODUCT DESCRIPTION

R 1	Succinate buffer ≥ 20 mmol/L; Bromocresol green ≥ 100 µmol/L; Stabilizer; detergent; Preservative.	X
	Phosphate Buffer ≥ 20 mmol/L; Bovine albumin at a concentration equivalent to 4,0 g/dL; Preservative. Traceable to reference material ERM-DA470k / IFCC	
STD		X

STABILITY IN USE

- The stability of the product (R1 and STD) in use is 24 months as long as followed by the recommended storage conditions (15 to 30 ° C).

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Reagent 1 (R1) and STD
Reagents are ready for use.

B) OPERATING RANGE

The operating range of the product is from 0,27 g/dL to 6,00 g/dL.

For higher values dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:

Calibrator Serum - Autocal H	REF	13.002.00
Normal Control Serum - Quantinorm		13.003.00
Pathological Controle Serum - QuantiAlt		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	5 µL	-
Sample	-	-	5 µL
R1	1,0 mL	1,0 mL	1,0 mL

- Homogenize well and keep the tubes for 10 minutes at a temperature of 15 to 30 ° C.
- Measure the Absorbance of the Standard and Sample versus White at 630 nm (600-640). The color is stable for 30 minutes.

B) CALCULATIONS

Albumin (g / dL) = $\frac{\text{Sample Absorbance} \times \text{Standard Concentration (g / dL)}}{\text{Standard Absorbance}}$

Example:

Standard Concentration = 4 g/dL

Sample Absorbance = 0.280

Standard Absorbance = 0.360

Albumin (g / dL) = $\frac{0,280 \times 4}{0,360} = 3.11$ g/dL

With Calibration Factor (CF):

CF = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Albumin (g / dL) = Sample Absorbance x Calibration Factor

Example:

Calibration Factor = $\frac{4}{0,360} = 11.1$

Albumin (g / dL) = 0.280 x 11.1 = 3.11 g/dL

Albumin (g / dL) = 0.280 x 11.1 = 3.11 g/dL

C) INTERPRETATION

Albumin is the most abundant plasma protein, synthesized primarily in the liver. Its main function is the maintenance of the colloidal osmotic pressure inside the blood vessels and extravascular spaces. It carries a large number of compounds: free fatty acids, phospholipids, metal ions, amino acids, drugs, hormones and bilirubin. Elevated albumin concentration occurs only in acute dehydration and is of no clinical importance. Its reduction occurs in several clinical conditions, such as: analbuminemia, inflammation, liver disease, urinary and gastrointestinal losses, energy and protein malnutrition, edema and ascites. Albumin can bind to numerous drugs, which is why low blood albumin concentrations have a significant pharmacokinetic effect

INTERFERING AND LIMITATIONS

Hemolyzed, Jaundice and Lipemic Sera: Hemoglobin > 500 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 2000 mg/dL interfere in the dosage.

Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 0,06 g/dL / quantification limit: 0,27 g/dL

Analytical Specificity: The product determines albumin specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 1,015x - 0,057$ and the correlation coefficient $r = 0,9979$. Using this equation the total systematic error estimated is -1,350% to a level of 2,0 g/dL and 0,075% to a level of 4 g/dL.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained.

Samples (g/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD (g/dL)	%CV	SD (g/dL)	%CV
2,2368	80	0,0153	0,7	0,0166	0,7
3,7708	80	0,0390	1,0	0,0505	1,3
5,0397	80	0,0284	0,6	0,0397	0,8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the material used are key factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Serum	3,5 – 5,2 g/dL
-------	----------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): g/L

Albumin (g/dL) x 10 = Albumin (g/L)

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer able to read at 630 nm (600 - 640).
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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