Alfa Amilase

Alpha Amylase | Alfa Amilasa
Ref. 11.001.00

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ENGLISH

INTENDED USE

Kit for the quantitative determination of the enzymatic activity of alpha-amylase in serum, plasma and urine. Use in vitro diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Alpha-amylase catalyzes the hydrolysis of α -(2-chloro-4-nitrophenyl)- β -1,4-galactopyranosylmaltoside (GalG2-CNP) into 2-chloro-4-nitrophenyl (CNP) and 1,4-galactopyranosylmaltoside (GalG2). The rate of CNP formation is proportional to the catalytic activity of the alpha-amylase in the sample.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum, plasma (heparin) and urine.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Amylase in serum and plasma is stable for 7 days if stored at 20 to 25 °C, 7 days if stored at 4 to 8 °C and one year stored at -20 °C. Urine amylase is stable for 2 days if stored at 20 to 25 °C, 10 days if stored at 4 to 8 °C and 3 weeks if stored at -20 °C.

PRODUCT DESCRIPTION

R 1

Buffer MES > 20 mM; potassium thiocyanate > 100mM; GalG2-CNP > 2 mM; activators; detergent; preservative.

X

Traceable to the specific absorptivity of the CNP chromogen using calibrated pipettes and handheld spectrophotometer that provide absolute values.

STABILITY IN USE

- After opening, the product (R1) in use is stable for 30 days, since recommended storage conditions (2 to 8 °C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

Ready-to-use reagent.

B) OPERATING INTERVAL

Serum / Plasma: The product's operating range is 4.2 U / L at 1200 U / L.

Urine: The operating range of the product is 5.2 U / L at 2000 U / L.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For calibration and internal control of laboratory quality it is recommended to use the serum calibrator and control sera below:

Calibrator Serum - Autocal H	13.002.00
Normal Control Serum - Quantinorm	13.003.00
Pathological Control Serum - Quantialt	13.004.00

REF

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette into test tube:

	Soro/Plasma	Urine
R1	1,0 mL	1,0 mL
Sample	20 µL	10 µL

2. Homogenize and insert into the cuvette holder thermostated at 37 °C, zeroed with purified water at 405nm. Set the stopwatch.

3. After 1 minute, note the initial absorbance A₀ and perform new readings every minute for 3 minutes (A₁, A₂ and A₃ respectively).

B) CALCULATIONS

Using the absorbance readings, calculate the average absorbance change per minute (ΔA / min):

$$(\Delta A/\text{min}) = \frac{(A_1 - A_0) + (A_2 - A_1) + (A_3 - A_2)}{3}$$

The amylase activity in the sample is calculated by multiplying the ΔA / min by the following factor:

Amilase (U/L)	
Soro/Plasma	$\Delta A/\text{min} \times 3292$
Urine	$\Delta A/\text{min} \times 6520$

Exemplo:

A₀ = 0,121 A₁ = 0,142

A₂ = 0,165 A₃ = 0,187

$$\Delta A/\text{min} = \frac{(0,142 - 0,121) + (0,165 - 0,142) + (0,187 - 0,165)}{3}$$

$\Delta A/\text{min} = 0,022$

Amilase (U/L) = 0,022 x 3292 = 72,4 U/L

C) INTERPRETATION

Alpha-amylase is a hydrolases-class enzyme that catalyzes the hydrolysis of 1,4- α -glycosidic bonds in polysaccharides. The activity of amylase present in serum and urine has both pancreatic and salivary gland origins, increasing greatly in acute pancreatitis and inflammation of this gland. Increased values are also found in several acute intra-abdominal disorders, in various extrahepatic conditions, and in neoplastic diseases.

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PERFORMANCE CHARACTERISTICS

Sensitivity:

Serum / Plasma: Limit of detection: 2.1 U / L / Limit of quantification: 4.2 U / L

Urine: Limit of detection: 1.26 U / L / Limit of quantification: 5.2 U / L

Analytical Specificity: The product specifically determines alpha-amylase in the presence of other interfering substances in the sample, up to the concentrations reported above..

Accuracy:

Serum / Plasma: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation Y = 0.992X - 1.702 and correlation coefficient R = 0.9916 was obtained. Using this equation the estimated total systematic error of -5.05% for a level of 40 U / L and -2.50% for a level of 100 U / L.

Urine: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation Y = 0.989X - 4.316 and correlation coefficient R = 0.9980 was

obtained. Using this equation the estimated total systematic error of -5.42% for a level of 100 U/L and -2.54% for a level of 300 U/L.

Precision:

Serum / Plasma: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (U/L)	Repetitions	In-run Precision		Total Precision	
		SD (U/L)	%CV	SD (U/L)	%CV
114,6	80	0,555	0,5	1,303	1,1
214,9	80	1,209	0,6	2,690	1,3
495,2	80	0,971	0,2	2,987	0,6

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

Urine: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained:

Samples (U/L)	Repetitions	In-run Precision		Total Precision	
		SD (U/L)	%CV	SD (U/L)	%CV
217,7	80	2,181	1,0	2,856	1,3
430,6	80	4,133	1,0	5,507	1,3
832,7	80	7,233	0,9	8,628	1,0

% CV: Coefficient of variation expressed as a percentage; SD: Standard deviation.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the use material are fundamental factors for the stability of the reagents and obtaining correct results.
- Additional care should be taken to avoid R1 contact with saliva or sweat.

REFERENCE RANGES

	Amylase (U/L)
Serum / Plasma	< 100
Urine	< 491

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): $\mu\text{Kat} / \text{L}$

$\text{Amylase (U/L)} \times 0,017 = \text{Amylase } (\mu\text{Kat/L})$

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 405 nm.
- Water bath, thermostated at 37 °C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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