Bilirrubina Direta e Total

Responsável Técnico:
Direct & Total Bilirubin| Dr. Gilson Sérgio Pizzo
Bilirrubina Directa y Total CRF MG – 5310
Ref. 10.003.00 **MS 80027310205**

FINALIDADE

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Recém-nascido	Bilirrubina Total (mg/dL)	
Idade	Prematuro	A termo
0-1 dia	1,0 a 8,0	2,0 a 6,0
1-2 dias	6,0 a 12,0	6,0 a 10,0
3-5 dias	10,0 a 14,0	4,0 a 8,0

	Adultos	mg/dL
Bilirrubina Direta	0,0 a 0,2	
Bilirrubina Total	0,0 a 2,0	

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ENGLISH

INTENDED USE

Kit for determination of Direct and Total Bilirubin in serum. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The determination is based on the diazo reaction, which colorimetrically quantifies the red staining of azobilirubin when bilirubin reacts under certain conditions with diazotized sulfanilic acid. The intensity of the formed color is proportional to the concentration of bilirubin in the sample. In the presence of dimethylsulfoxide (accelerator), the indirect form (unconjugated) and the direct form (monoconjugate, diconjugate and protein bound) are reacted to determine Total Bilirubin. The direct form is determined without the accelerator.

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Immediately after collection, protect the sample from the action of light.

Preservation: Serum bilirubin is stable for 7 days at 4 to 8 °C and 6 months at -20 °C.

PRODUCT DESCRIPTION

R BD	Sulfamic acid ≥ 10mmol / L, hydrochloric acid ≥ 0.5% v / v.	
R BT	Sulfamic acid ≥ 10mmol / L, hydrochloric acid ≥ 0.1% v / v, Dimethylsulfoxide ≥ 30% v / v.	
R NI	Sodium nitrite ≥ 10 mmol / L, sucrose ≥ 2 mol / L, preservative.	

STABILITY IN USE

- Once opened, the product (RBD, RBT and RNI) in use is stable to the expiration date printed on the label, since recommended storage conditions (2 to 8 °C) are followed.
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

RBD, RBT e RNI: Ready-to-use reagents.

B) OPERATING INTERVAL

The operating range of Direct Bilirubin is from 0.17 mg / dL to 15.0 mg / dL.

The operational range of Total Bilirubin is from 0.18 mg / dL to 15.0 mg / dL.

For higher values, dilute the sample with 150 mL NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Laboratory Quality Control it is recommended to use the calibrator serum and the control sera below:

Serum Calibrator - Autocal H	13.002.00
Normal Control Serum- Quantinorm	REF 13.003.00
Serum Pathological Control - Quantialt	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) DIRECT BILIRRUBIN

1 - Pipette into test tube.

	Blank	Blank Sample	Sample
RBD	1,0 mL	1,0 mL	1,0 mL
RNI	---	---	30 µL
Sample/Calibrator	---	50 µL	50 µL

2. Homogenize well and keep the tubes for 5 minutes at room temperature.

3. Measure the absorbance of the Calibrator, Calibrator Blank, Sample and Sample Blank against white at 550 nm (530-570). Color is stable for 5 minutes.

B) TOTAL BILIRRUBIN

1. Pipette into test tubes:

	Blank	Blank Sample	Sample
RBT	1,0 mL	1,0 mL	1,0 mL
RNI	---	---	30 µL
Sample/Calibrator	---	50 µL	50 µL

2. Homogenize well and keep the tubes for 5 minutes at room temperature.

3. Measure the absorbance of the Calibrator, Calibrator White, Sample and Sample Blank against white at 550 nm (530-570). Color is stable for 30 minutes.

C) CALIBRATION

Use Calibrator Autocal H Ref.13.002.00. The concentration of Bilirubin in the calibrator is traceable to the Standard Reference Material (SRM) 916a of the National Institute of Standards and Technology (NIST).

D) CALCULATIONS

Bilirubin (mg / dL) = $\frac{(\text{Sample Abs} - \text{Sample White Abs})}{(\text{Calibrator Abs} - \text{White Calibrator Abs})}$ x Calibrator Value (mg / dL)

Exemple:

Calibrator Concentration = 4,62 mg / dL

Sample Absorbance = 0,194

Sample Blank Absorbance = 0,040

Calibrator Absorbance = 0,234

Absorbance Blank Calibrator = 0,024

Bilirubin (mg / dL) = $\frac{0,194 - 0,040}{0,234 - 0,024}$ x 4,62 = 3,39 mg / dL

With Calibration Factor

Calibration Factor = $\frac{\text{Calibrator Concentration (mg / dL)}}{(\text{Calibrator Abs} - \text{White Calibrator Abs})}$

Bilirubin (mg / dL) = (Sample Abs - Abs Sample White) x Calibration Factor

Exemple:

Calibration Factor = $\frac{4,62}{0,234 - 0,024}$ = 22

Bilirubin (mg / dL) = (0,194-0,040) x 22 = 3,39 mg / dL

C) INTERPRETATION

Bilirubin is a pigment derived from heme degradation. Bilirubin is transported to the liver, poorly associated with albumin, in its unconjugated form. Carried into the hepatocyte it is rapidly conjugated to produce the bilirubin glucuronides, which are then excreted into the bile. Elevations of conjugated (direct) bilirubin are usually specific markers of hepatic dysfunction, acute hepatitis, or cholestasis. Increased unconjugated (indirect) bilirubin is usually caused by hemolytic processes or by a reduction in liver metabolism.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin \geq 150 mg / dL, hemolysis interferes in the dosages of direct bilirubin and total bilirubin.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity:

	Direct Bilirubin	Total Bilirubin
Detection Limit	0,14 mg/dL	0,13 mg/dL
Qualification Limit	0,17 mg/dL	0,18 mg/dL

Analytical Specificity: The product specifically determines Direct and Total Bilirubin in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The results are as follows:

Direct Bilirubin: Regression equation $y = 1,017x - 0,008$ and correlation coefficient $r = 0,9909$. Using this equation the estimated total systematic error of 0,9% for a level of 1,0 mg / dL and 1,6% for a level of 8,0 mg / dL.

Total Bilirubin: Regression equation $y = 1,014x - 0,010$ and correlation coefficient $r = 0,9901$. Using this equation the total systematic error estimated was 1,20% for a level of 5,0 mg / dL and 1,32% for a level of 13,0 mg / dL.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (mg/dL)	Repetitions	Within-run precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
0,767	80	0,032	4,2	0,039	5,1
3,455	80	0,072	2,1	0,181	5,2
12,218	80	0,067	0,5	0,172	1,4

Samples (mg/dL)	Repetitions	Within-run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
1,0471	80	0,0135	1,3	0,0135	1,3
5,4074	80	0,0519	1,0	0,0623	1,2
13,4034	80	0,0463	0,3	0,1097	0,8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Newborn		Total Bilirubin (mg/dL)	
Age	Premature	Normal	
0-1 day	1,0 to 8,0	2,0 to 6,0	
1-2 day	6,0 to 12,0	6,0 to 10,0	
3-5 day	10,0 to 14,0	4,0 to 8,0	

Adults		mg/dL
Direct Bilirubin		0,0 to 0,2
Total Bilirubin		0,0 to 2,0

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 550 nm (530-570).
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.bioteecnica.ind.br or by telephone + 55 35 3214 4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use Biotécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. Adaptations of semi-automatic protocols with preparation of Working Reagents, data concerning the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory Board of Biotécnica Ltda, by phone +55 35 3214 4646 or by email sac@biotecnicatda.com.br.

AUTOMATION

This procedure is automated in most analyzers. The applications are available at www.biotecnicatda.ind.br

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