## Bilirrubina Direta

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### ENGLISH

#### INTENDED USE

Kit intended to determine direct bilirubin in serum and plasma (heparin and EDTA) samples. Diagnostic use only.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C and protect from light. The product must remain out of the specified temperature only the time required for testing.
- Reagents ready for use.
- Once opened, the product is stable until the expiration date printed on the label, as long as the recommended storage conditions (2 to 8°C) are followed.
- Do not use reagents whose shelf life has expired.

#### WORKING PRINCIPLE

**Method:** DCA (Diazotized Dichloroaniline).

At acidic pH, direct bilirubin is coupled with diazotized dichloroaniline to form a colored complex, known as azobilirubin, which can be measured at 550 nm (500-590 nm). The color intensity is proportional to the concentration of direct bilirubin present in the sample.

#### SAMPLE: TYPE, COLLECTION, HANDLING AND STABILITY

**Sample Type:** serum and plasma (heparin and EDTA).

**Collection and Handling:** collect the sample in accordance with the Good Laboratory Practices. All samples should be treated as potentially infectious material.

#### Preservation:

	Temperature	Stability Period
Serum and Plasma*	4 to 8 °C	7 months
	-20 °C	6 months

\*Necessary to keep out of the light when stored for longer than 8 hours.

### PRODUCT DESCRIPTION

<b>R 1</b>	Sulfamic Acid <300 mmol/L and stabilizers.
<b>R 2</b>	Diazotized dichloroaniline <2.0 mmol/L and stabilizer. pH <5

### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. For the internal laboratory quality control, it is recommended the use of the calibrator and controls below:

Calibrador - Autocal H	13.002.00
Normal Control - Qualitorm	13.003.00
Pathological Control - Quantialt	13.004.00

**REF**

### NECESSARY EQUIPMENT FOR TESTING

- Spectrophotometer or photometer for reading at 550 nm (505 – 590 nm).
- Thermostatic water bath at 37 °C and test tubes.
- Glass pipettes and/or automatic, clock or chronometer.

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Pipette in the test tubes:

	Blank	Sample Blank	Sample	Calibrator Blank	Calibrator
Calibrador	---	---	---	80 µL	80 µL
Sample	---	80 µL	80 µL	---	---
Purified water	80 µL	---	---	---	---
R1	1,0 mL	1,0 mL	1,0 mL	1,0 mL	1,0 mL
R2	250 µL	---	250 µL	---	250 µL

Homogenize and incubate the test tubes for 1 minute at 37°C

R2	250 µL	---	250 µL	---	250 µL
----	--------	-----	--------	-----	--------

2. Homogenize and incubate at 37 °C for 7 minutes.

3. Measure the Blank's, Sample Blank's, Sample's and Standard's absorbance at 550 nm (500-590 nm). The final reaction is stable for 30 minutes.

### B) CALCULATIONS

Direct Bilirubin (mg/dL) =  $\frac{\text{Sample's } \Delta A - \text{Blank}}{\text{Calibrator's } \Delta A - \text{Blank}}$  x Calibrator Concentration (mg/dL)

In which:

$\Delta A$  = Sample's or Calibrator's Absorbance – Sample or Calibrator Blank's Absorbance

#### Calculations with the Calibration Factor:

Calibration Factor =  $\frac{\text{Calibrator Concentration (mg/dL)}}{\text{Calibrator's } \Delta A - \text{Blank}}$

Direct Bilirubin (mg/dL) =  $(\text{Sample's } \Delta A - \text{Blank}) \times \text{Calibration Factor}$

**Automation:** this product is compatible to most types of automatic analyzers. Instrument settings are available at [www.biotecnicaltda.ind.br](http://www.biotecnicaltda.ind.br)

### C) INTERPRETATION

Bilirubin is conjugated in the liver and excreted in bile, resulting from the breakdown of hemoglobin in the reticuloendothelial system. Its dosage assists in the diagnosis of hepatobiliary and hemolytic diseases. There are inherited or acquired diseases that are responsible for altering the stages involved in its production, uptake, storage, metabolism and excretion, changing the balance between unconjugated (indirect), conjugated (direct) or both (total) bilirubin, causing hyperbilirubinemia. Direct bilirubin increases in the body due to inherited liver diseases such as Dubin-Johnson and Rotor diseases, damage to hepatocytes (viral, toxic or alcoholic), biliary obstruction (lithiasis or neoplasms), acute or chronic hepatitis and toxic reactions to various medications (such as chlorpromazine).

### PERFORMANCE CHARACTERISTICS

Operating range
0.06 to 20.00 mg/dL

For concentrations above the operating range, dilute the sample with NaCl 150 mM (0.9%), proceed with a new dosage and multiply the result by the dilution factor.

Sensitivity	
Detection Limit	Quantification Limit
0.04 mg/dL	0.06 mg/dL

Analytical Specificity		
Hemoglobin	Triglycerides	Ascorbic Acid
200 mg/dL	750 mg/dL	8 mg/dL

Interfering substances up to the values presented above do not cause significant alterations in results. For drugs, consult the recommended reference (Young, 2000).

Accuracy	
Number of Samples	40 in duplicate
Regression Equation	y = 0.9228x + 0.0565
Correlation Coefficient (R)	0.9992

**Serum:** by applying the regression equation, the total systematic error estimated for decision levels of 1.40 mg/dL, 2.50 mg/dL and 20.00 mg/dL was -3.7%, -5.5% and -7.4%, respectively.

#### Precision:

Within-run precision was determined with 5 replicates for 10 days; between-run precision was determined in 3 condition, with 3 replicates for 5 days.

Samples (mg/dL)	Within-Run Precision SD (mg/dL) CV (%)	Samples (mg/dL)	Between-Run Precision SD (mg/dL) CV (%)
0.60	0.02 4.00	0.57	0.02 3.80
4.18	0.10 2.50	4.17	0.07 1.60
7.25	0.11 1.50	7.25	0.10 1.40

CV: Coefficient of variation; SD: Standard deviation.

### RESIDUAL RISKS, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with the Good Laboratory Practices.
- Follow the Good Laboratory Practices' instructions to establish the quality of water.
- Do not mix reagents from different lots or exchange the caps from different reagents in order to avoid cross contamination. Do not use the reagent if it displays any signs in disagreement with the ones specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions.
- The level of the water bath must be greater than that of the test tubes containing the reaction.

### REFERENCE RANGES

All ages	Up to 0.4 mg/dL
----------	-----------------

These values are intended for orientation only. It is recommended that each laboratory establishes its own reference ranges.

## WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to the Good Laboratory Practices, in a proper place for potentially infectious material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or calling for +55 35 3214 4646

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

All Biotécnica products are made according to the Good Manufacturing Practices and other current sanitary regulations. Their performance is assured as long as all Biotécnica instructions are followed. In case of doubt while using the product, contact our Scientific Advisory team by calling +55 35 3214 4646, your local distributor or sending an e-mail to [sac@biotechnicaltda.com.br](mailto:sac@biotechnicaltda.com.br).

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