Cálculo

Calcium | Calcio
Ref. 12.001.00

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ENGLISH

INTENDED USE

Kit intended to determination of calcium in serum, plasma and urine. Diagnostic use only.

STORAGE AND HANDLING

- Store at 15 to 30 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Calcium forms with ortho-cresolphthalein complexone, in alkaline medium, a colored complex that can be measured at 550-590 nm. The intensity of the formed color is proportional to the amount of calcium present in the sample.

SAMPLE – PREPARATION AND STABILITY

Sample Type: Serum, plasma (heparin) and urine.

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The calcium in serum and plasma (heparin) is stable for 14 days if stored at a temperature of 4 to 8 °C and 6 months if stored at a temperature of -20 °C. Urine: Use sample collected within 24 hours in a bottle containing 20 to 30 mL of 6 mol / L HCl solution to prevent the precipitation of calcium salts. Mix well, measure the volume and proceed with the test.

PRODUCT DESCRIPTION

| | | | |
|------------|-----------------------------------------------------------------------------------------------------------------------------------|--|--|
| R 1 | 2-Amino methyl propanol> 1% v/v; Monoethanolamine> 1% v/v and stabilizers. | | |
| R 2 | Hydrochloric acid> 0.1% v / v; Ortho-cresolphthalein complexone> 0.001% w / v; 8-hydroxyquinoline > 1% mmol / L and preservative. | | |
| STD | Calcium chloride in concentration equivalent to 10 mg/dL and preservative. Traceable to reference material NIST 3109a. | | |

STABILITY IN USE

- The stability of the product (R1,R2 and STD) in use is 24 months as long as followed by the recommended storage conditions (15 to 30 ° C).

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

Working Reagent:

Mix in proportion: 1 part of R1 + 1 part of R2. Homogenize gently. Stable for 24 hours at 15 to 30 ° C.

NOTE: The absorbance of the working reagent should be less than 0.200 at 570 nm in equipment cleared with purified water.

STD: Ready for use.

B) OPERATING INTERVAL

The operating range of the product is from 0.91 mg / dL to 18.00 mg / dL.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor...

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:

| | |
|-----------------------------------------|-----------|
| Calibrator Serum - Autocal H | 13.002.00 |
| Normal Control Serum – Quantinorm | 13.003.00 |
| Pathological Controle Serum – QuantiAlt | 13.004.00 |

REF

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the assay tubes:

| | Blank | Standard | Sample |
|--------|--------|----------|--------|
| STD | - | 10 µL | - |
| Sample | - | - | 10 µL |
| R1 | 1,0 mL | 1,0 mL | 1,0 mL |

2. Homogenize and incubate 2 minutes at 15 to 30 °C.

3. Measure the absorbance of Standard (Ast) and Sample (Asa) against the blank at 570 nm (550-590 nm). The final reaction is stable for 30 min.

B) CALCULATIONS

Calculation of serum

Calcium (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{STD Concentration}$

Example:

Standard Concentration = 10 mg / dL

Sample Absorbance = 0,220

Standard Absorbance = 0,250

Calcium (mg / dL) = $\frac{0,220}{0,250} \times 10 = 8,8$ mg / dL

With Calibration Factor (CF):

CF = $\frac{\text{STD Concentration (mg/dL)}}{\text{STD Absorbance}}$

Calcium (mg/dL) = Sample Absorbance x CF

Example:

Calibration Factor = $\frac{10}{0,250} = 40$

Cálcium (mg/dL) = 0,220 x 40 = 8,8 mg/dL

Urine 24 hours:

Calcium urine (mg / dL) = $\frac{\text{Sample Absorbance}}{\text{Standard Absorbance}} \times \text{Standard Concentration (mg / dL)}$

Calcium (mg/24h) = $\frac{\text{urine calcium (mg/dL)} \times \text{urinary volume of 24 hr (mL)}}{100}$

C) INTERPRETATION

Calcium is the most predominant cation in the body. The skeleton contains 99% of body calcium with an unknown structure and composition roughly that of hydroxyapatite. The remainder of the calcium is distributed in tissues and extracellular fluid where it performs specific functions in vital processes, among them: blood coagulation, neuromuscular conduction, cardiac and skeletal muscle excitability, enzymatic activation, and the preservation of cell membrane permeability and integrity. Serum calcium levels are controlled by parathyroid hormone (PTH), calcitonin, and vitamin D. An imbalance in some of these modulators causes changes in serum and body calcium levels. Serum increases in PTH or vitamin D are generally associated with hypercalcemia. Increased levels of serum calcium can be observed in multiple myeloma and other neoplastic diseases. Hypocalcemia can be observed for example in hypoparathyroidism, nephrosis and pancreatitis.

INTERFERING AND LIMITATIONS

Hemolyzed, Jaundice and Lipemic Sera: Hemoglobin > 500 mg/dL, Bilirubin > 35 mg/dL, Triglyceride > 1000 mg/dL. Interfere in the dosage

Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 0,58 mg/dL / quantification limit: 0,91 mg/dL

Analytical Specificity: The product determines calcium specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,990x + 0,090$ and the correlation coefficient $r = 0,9984$. Using this equation the total systematic error estimated is 0,80% to a level of 5,0 mg/dL and -0,25% to a level of 12,0 mg/dL.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained:

| Samples (mg/dL) | Repetitions | Within-Run Precision | | Total Precision | |
|-----------------|-------------|----------------------|-----|-----------------|-----|
| | | SD | %CV | SD | %CV |
| 10,2688 | 80 | 0,0564 | 0,5 | 0,1332 | 1,3 |
| 13,2685 | 80 | 0,0229 | 0,2 | 0,1021 | 2,1 |
| 16,9265 | 80 | 0,0118 | 0,1 | 0,1542 | 0,9 |

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

| | mg/dL | mmol/L |
|------------------|-------------|---------------|
| Serum and Plasma | 8,6 - 10,3 | 2,15 - 2,57 |
| | mg/24 hours | mmol/24 hours |
| Urine | Up to 300 | Up to 7,49 |

These values are intended for orientation only. It's recommended that each lab establishes its own reference range. Conversion for International System of Units (SI): mmol/L Calcium (mg/dL) X 0,25 = Calcium (mmol/L)

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer able to read at 570 nm (550-590).
- Pipettes and micropipettes.
- Clock or chronometer.
- Assay tubes.

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for +55 35 3214 4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels

of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotecnicaltda.ind.br

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