Desidrogenase Láctica

Lactate Dehydrogenase (LDH) | Lactato Deshidrogenasa

Ref. 11.004.00

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ENGLISH

INTENDED USE


Kit intended to determination of Lactate Dehydrogenase activity in serum and plasma. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The Lactate Dehydrogenase (LDH) catalyzes the reduction of pyruvate to lactate in the presence of NADH. The rate of conversion of NADH to NAD⁺ in the reaction is proportional to the catalytic activity of LDH determined by the measurement of the absorbance decrease at 340 nm.

SAMPLE – PREPARATION AND STABILITY

Sample Type: Serum and plasma (heparin)

Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Separate the serum / plasma from the cell elements immediately. Do not use hemolyzed sample.

Preservation: The Lactate Dehydrogenase in serum and plasma (heparin) is stable for 1 day if stored at a temperature of 4 to 8 °C and 3 weeks at a temperature of -20 °C.

PRODUCT DESCRIPTION

R 1	TRIS buffer ≥ 50 mmol/L; Sodium pyruvate ≥ 0,5 mmol/L; activators; stabilizers and preservative.
R 2	Carbonate buffer ≥ 2 mmol/L; NADH ≥ 0,5 mM; preservative.

Optimized and standardized method according to the recommendations of the Deutschen Gesellschaft für Klinische Chemie (DGKC), using calibrated pipettes, manual spectrophotometer and specific absorptivity of the chromogen that provide absolute values.

STABILITY IN USE

- The stability of the product (R1 and R2) in use is 18 months as long as followed by the recommended storage conditions (2 to 8 °C).
- The stability of working reagent is 30 days, as long as followed by the conditions of preparation and recommended storage (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE

A) REAGENT PREPARATION

R1 and R2: Ready-to-use reagents.

Working Reagent: Prepare the Working Reagent in the proportion of: 4 parts of R1 + 1 part of R2. Homogenize gently. Stable for 30 days at 2 to 8 °C.

B) OPERATING RANGE

The product operating range is from 11,38 U/L to 2000 U/L. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:

Calibrator Serum - Autocal H	REF	13.002.00
Normal Control Serum – Quantinorm		13.003.00
Pathological Controle Serum – QuantiAlt		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Heat the reagent for 3 minutes at 37 °C.
2. Pipette in the assay tube:

Work Reagent	Volume
Sample	1.0 mL 20 µL

3. Homogenize and insert it immediately in the thermostated cuvette at 37 °C.
4. After 1 min. make a note of the initial absorbance (A₀) and read again after exactly 1, 2 and 3 minutes (A₁, A₂ and A₃ respectively).

B) CALCULATIONS

Using the measured absorbances, calculate the mean variation of absorbance per minute (ΔA/min):

$$\Delta A/\text{min.} = \frac{(A_0 - A_1) + (A_1 - A_2) + (A_2 - A_3)}{3}$$

The activity of the LDH is then calculate by multiplying the ΔA/min by the factor below:

$$\text{LDH (U/L)} = \Delta A/\text{min} \times 8095$$

Example:

A₀ = 1.310 - A₁ = 1.265 - A₂ = 1.220 - A₃ = 1.175

$$\Delta A/\text{min} = \frac{(1,310 - 1,265) + (1,265 - 1,220) + (1,220 - 1,175)}{3}$$

ΔA/min = 0.045

LDH (U/L) = 0.045 x 8095 = 364 U/L

C) INTERPRETATION

LDH has five isoenzymes, which can be separated according to their electrophoretic mobility. Each isoenzyme is a tetramer composed of two subunits M and / or H. It is widely distributed in tissues, particularly heart, liver, muscle and kidneys at concentrations approximately 500 times greater than those normally found in serum. Therefore, it can be used as a marker of acute or chronic tissue injury and, sometimes, in the follow-up of progressive lesions. Serum LDH elevation occurs in several clinical conditions, such as myocardial infarction, hemolysis, neoplasms, liver, kidney, lung and muscle disorders.

INTERFERING AND LIMITATIONS

Hemolyzed, Jaundice and Lipemic Sera: Hemoglobin > 25 mg/dL, Bilirubin > 40 mg/dL, Triglyceride > 1500 mg/dL.

Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 5,22 U/L / quantification limit: 11,38 U/L.

Analytical Specificity: The product determines Lactate Dehydrogenase specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.

Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,999x - 6,54$ and the correlation coefficient $r = 0,9958$. Using this equation the total systematic error estimated is -1,73% to a level of 400 U/L and -0,64% to a level of 1200 U/L.

Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained.

Samples (U/L)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
261,059	80	2,709	1,0	2,834	1,1
776,150	80	2,239	0,3	2,496	0,3
1320,326	80	8,64	0,7	9,191	0,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The level of water in the thermostatic bath must be over the level of mix reaction within the assay tubes.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES

ADULTS	37 °C
Men	214 - 450 U/L
Women	195 - 453 U/L

CHILDREN	37 °C
0 a 6 months	613 - 1020 U/L

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI):

$$LDH (U/L) \times 0,017 = LDH (\mu\text{kat})$$

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 340 nm.
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotechnica.ind.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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