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#### ENGLISH

**INTENDED USE**  
Kit for determination of phosphorus in serum and urine. Use in vitro diagnostic.

#### STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Keep out of the light.
- Stable up to the expiration date of the kit that is printed on the package label.
- Do not use reagents whose expiration date has expired.

#### WORKING PRINCIPLE

The inorganic phosphate in the sample reacts with the ammonium molybdate in an acid medium, giving rise to the phosphomolybdate complex, which can be measured at 340 nm. The color intensity is proportional to the concentration of inorganic phosphate in the sample.

Inorganic phosphate + ammonium molybdate  $\xrightarrow{H_2SO_4}$  Phospho-molybdenum complex

#### SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

**Sample Type:** Serum and Urine.  
**Collection, handling and preparation:** Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material. Serum samples should be separated from the cell elements immediately.  
**Preservation:** in the serum phosphorus is stable for 1 day if stored at 20-25 °C, 4 days at 4 to 8 °C and 1 year at -20 °C. In urine the phosphorus is stable for 2 days if stored at 20 to 25 °C and 6 months at 4 to 8 °C at pH <5.0.  
**Urine:** Use sample collected within 24 hours in a vial containing 20 to 30 mL of HCl 50% to prevent phosphate precipitation. Mix well, measure the volume, check the pH and if necessary adjust it to below 3.0. Centrifuge samples containing precipitates and use the supernatant. Dilute an aliquot of urine in the proportion of 1:10 with purified water (where necessary, dilution of the urine should be altered to achieve a result within the operating range of the method). Use diluted urine for testing. Multiply the result obtained by 10 (or by the factor that has been applied).

#### PRODUCT DESCRIPTION

<b>R 1</b>	Ammonium molybdate ≥ 0.5 mmol/L; Sulfuric acid ≥ 0.5% v/v; Triton X100 ≥ 0.1% v/v.	
<b>STD</b>	Sodium phosphate in concentration equivalent to 5.0 mg/dL, preservative. Traceable to reference material NIST 3139a.	

#### STABILITY IN USE

- The stability of the product (R1 and STD) in use is 24 months, since the recommended storage conditions (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

#### TREATMENT AND HANDLING OF THE PRODUCT TECHNICAL PROCEDURE

**A) REAGENT PREPARATION**  
R1 and STD: Ready for use.

**B) OPERATING INTERVAL**  
**Serum:** 0.35 mg dL to 15.00 mg/dL.  
For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.  
**Urina:** 6.86 mg / dL at 233.85 mg / dL, corrected for dilution factor

#### QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and the control sera below:

Serum Calibrator - Autocal H	13.002.00
Normal Control Serum - Quantinorm	13.003.00
Serum Pathological Control - Quantial	13.004.00

**REF**

#### TEST PROCEDURE, CALCULATION AND INTERPRETATION

**A) TEST PROCEDURE**  
1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 5 minutes at 37 °C.
3. Measure the absorbance of the Standard and the Sample against the reagent Blank at 340 nm. The final color of the reaction is stable for 60 minutes.

#### B) CALCULATIONS

**Serum:**  
Phosphorus (mg/dL) =  $\frac{\text{Sample Absorbance} \times \text{Standard Concentration (mg/dL)}}{\text{Standard Absorbance}}$

**Exemplo:**  
Standard Concentration = 5.0 mg / dL  
Sample Absorbance = 0.225  
Standard Absorbance = 0.252  
Phosphorus (mg/dL) =  $\frac{0.225 \times 5.0}{0.252} = 4.46$  mg / dL

**With Calibration Factor:**

Calibration Factor =  $\frac{\text{Standard Concentration (mg/dL)}}{\text{Standard Absorbance}}$

Phosphorus (mg/dL) = Sample Absorbance x Calibration Factor

**Example:**

Calibration Factor =  $\frac{5,0}{0,252} = 19,8$

Phosphorus (mg / dL) = 0.225 x 19.8 = 4.46 mg/dL

**24-hour urine:**

Phosphorus Urine (mg/dL) =  $\frac{\text{Sample Absorbance} \times \text{Standard Concentration (mg / dL)} \times 10^*}{\text{Standard Absorbance}}$

\* 10 = Dilution factor applied.

Phosphorus urine (mg/24 hours) =  $\frac{\text{Phosphorus urine (mg/dL)} \times 24\text{-hour urine volume (in mL)}}{100}$

**C) INTERPRETATION**

Inorganic phosphate is located in the bone in the form of hydroxyapatite. In the extracellular fluid, phosphorus occurs as inorganic phosphate and in small amounts as organic phosphate in the form of phospholipids. Organic phosphate is found in cells such as phospholipids, nucleic acids and adenosine triphosphate - ATP. The mechanism of regulation of the levels of phosphorus and calcium ions in the blood are influenced by the interactions between parathyroid hormone and vitamin D. Phosphorus is increased in serum in cases of renal failure, hypoparathyroidism, resistance to parathyroid hormone and acromegaly. Its decrease occurs in vitamin D deficiency, hyperparathyroidism and Fanconi's syndrome

**INTERFERING AND LIMITATIONS****Hemolysis, Jaundice and Lipemia:**

Serum: Hemoglobin > 200 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 350 mg/dL.

Urine: Hemoglobin > 180 mg/dL / Bilirubin > 20 mg/dL.

Medicines: consult recommended literature reference (Young, 2000).

**PERFORMANCE CHARACTERISTICS****Sensitivity:**

Serum: Detection limit: 0.08 mg/dL / Limit of quantification: 0.35 mg/dL.

Urine: Limit of detection: 0.050 mg/dL / Limit of quantification: 6.86 mg/dL, corrected by dilution factor.

**Analytical Specificity:** The product specifically determines phosphorus in the presence of other interfering substances in the sample, up to the concentrations reported above.

**Accuracy:**

Serum: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation  $y = 0,989x + 0,069$  and correlation coefficient  $r = 0,9943$  was obtained. Using this equation the estimated total systematic error of 1.20% for a level of 3.00 mg/dL and 0.28% for a level of 5.00 mg/dL

Urine: The method was compared with a similar method by the determination of 20 samples in duplicate. The regression equation  $y = 1,006x - 2,703$  and correlation coefficient  $r = 0,9995$  was obtained. Using this equation, the total systematic error estimated was -0.17% for a level of 350.0 mg/24hs and 0.29% for a level of 900.0 mg/24hs.

**Precision:**

Serum: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (mg/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD (mg/dL)	%CV
3,64	80	0,0414	1,1	0,0522	1,4
7,61	80	0,0063	0,1	0,0924	1,2
14,53	80	0,0793	0,5	0,1266	0,9

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

Urine: It was determined using samples at 02 decision levels, two runs in duplicate per day for 20 days, being obtained:

Samples (mg/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD (mg/dL)	%CV	SD	%CV
26,93	80	0,18	0,7	0,22	0,8
87,76	80	0,38	0,4	0,58	0,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

**RESIDUAL RISK, WARNINGS AND PRECAUTIONS**

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

**REFERENCE RANGES**

	Serum	Urine
Adults	2,5 - 4,5 mg/dL	400 - 1300 mg/24 horas
Children	4,0 – 7,0 mg/dL	

These values are intended for orientation only. It's recommended that each lab establishes its own reference ranges.

Conversion for International System of Units (SI): mmol/L

Phosphorus (mg/dL) x 0,323 = Phosphorus (mmol/L)

**ADDITIONAL EQUIPMENT**

- Spectrophotometer or photometer able to read at 340 nm
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

**WARNINGS AND PRECAUTIONS**

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at [www.biotecnica.ind.br](http://www.biotecnica.ind.br) or calling for 55-35-3214-4646

**QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE**

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: [sac@biotecnica.com.br](mailto:sac@biotecnica.com.br)

**AUTOMATION**

This procedure is automated in most analyzers. The applications are available at [www.biotecnica.com.br](http://www.biotecnica.com.br)

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- YEE, H.Y. A simplified method for automated phosphorus analysis in serum and urine. *Clin. Chem.* v.14, p.898-903, 1968.
- DALY, A. J.; ERTINGSHAUSEN, G. Direct method for determining inorganic phosphate in serum with the "CentrifChem". *Clin. Chem.* v.18, p.263-265, 1972.
- GAMST, O.; TRY, K. Determination of serum phosphate without deprotonization by ultraviolet spectrophotometry of the phosphomolybdic acid complex. *Scand. J. Clin. Lab. Invest.* v.40, p.483-486, 1980.
- BURTIS, Carl A.; ASHWOOD, Edward R.; BRUNS, David E. *Tietz: Fundamentos de Química Clínica.* 6. ed. Rio de Janeiro: Saunders Elsevier, 2008. 959 p.
- YOUNG, D. S. Effects of drugs on clinical laboratory tests - vol. 2, 5 ed. Washington DC: AAC Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981.

**TABELA DE SÍMBOLOS INTERNACIONAIS / TABLE OF INTERNATIONAL SYMBOLS /****TABLA DE SÍMBOLOS INTERNACIONALES**

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