FERRO CRX

Iron CRX | Hierro CRX
Ref. 12.004.00

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ENGLISH

INTENDED USE

Kit intended for determination of iron in the serum. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 15 to 30 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

The iron present in the sample, at acidic pH, is released from transferrin. Fe³⁺ ions are reduced by hydroxylamine to Fe²⁺. These react with the chromogen ferrozine to form a colored complex, which exhibits a maximum absorption at 560 nm. The intensity of the color is proportional to the concentration of iron in the sample.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practices. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The serum iron is stable for 7 days if stored at a temperature of 20 to 25 °C, 3 weeks if stored at a temperature of 4 to 8 ° C and 1 year if stored at -20 ° C.

PRODUCT HANDLING

R 1 Acetate Buffer ≥ 60 mmol / L; Cetyltrimethylammonium bromide ≥ 0.1 mmol / L; Cromazurol B ≥ 50 μmol / L; preservative; stabilizers.

STD Solution of Fe3 + ions in concentration equivalent to 100 μg /dL, activator, reducer; stabilizer. Traceable to reference material NIST -3126a and NIST-928.

STABILITY IN USE

- After opening, the product (R1 and STD) in use is stable until the expiration date printed on the label, since recommended storage conditions are followed (15 to 30 ° C).

Reagents should remain outside the specified temperature only for the time required to perform the tests.

PRODUCT TREATMENT AND HANDLING

A) REAGENT PREPARATION

R1 and STD: Ready-to-use reagents.

B) OPERATING INTERVAL

The operating range of the product is from 11.74 μg/dL to 500.24 μg/dL.

For higher values, dilute the sample with 150 mmol/L NaCl (0.9% w/v), carry out a new dosage and multiply the result obtained by the dilution factor .

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Control of Laboratory Quality it is recommended to use the calibrator serum and the control sera below:

Serum Calibrator - Autocal H	13.002.00
Normal Control Serum – Quantinorm	REF 13.003.00
Serum Pathological Control– Quantialt	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the assay tubes:

	White	STD	Sample
Sample	----	----	50 μL
STD	----	50 μL	----
R1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize gently and incubate for 5 minutes at a temperature of 10 to 30 ° C.

3. Measure the absorbance of the Standard and the Sample against the White at 620 nm (620-630 nm). Color is stable for 45 minutes in the dark.

B) CALCULATIONS

Iron (μg/dL) = $\frac{\text{Sample absorbance}}{\text{Standard absorbance}} \times \text{Standard Concentration (μg/dL)}$

Example:

Standard Concentration = 200 μg/dL

Absorbance of the sample = 0,161

Absorbance of standard = 0,230

Iron (μg/dL) = $0,161 \times \frac{200}{0,230} = 140,0 \mu\text{g/dL}$

With calibration factor:

Calibration factor = $\frac{\text{Standard Concentration (μg/dL)}}{\text{Standard Absorbance}}$

Iron (μg/dL) = Sample Absorbance x Calibration factor

Example:

Calibration factor = $\frac{200}{0,230} = 869,56$

Iron (μg/dL) = 0,161 x 869,56 = 140,0 μg/dL

C) INTERPRETATION

In the duodenum and upper jejunum, absorption of iron ingested as Fe²⁺ occurs. The free trivalent form and Fe³⁺ bound to the hemeproteins present in the diet are reduced by vitamin C. About 1 mg of iron is absorbed daily. Once in the mucosal cells, Fe²⁺ ions bind to transporters. Before passing into plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bind to apotransferrin. The apotransferrin-Fe³⁺ complex is called transferrin which carries the iron ions into the blood plasma. The determination of iron (non-heme) is used in the diagnosis and treatment of diseases such as anemia due to iron deficiency, hemochromatosis, chronic kidney disease and in the monitoring of other types of anemias.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Bilirubin> 40 mg/dL lipemia e hemolysis interfere with the dosage.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 9,72 μg/dL / Limit of quantification: 11,74 μg/dL.

Analytical Specificity: The product specifically determines iron in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation y = 0,990x + 0,755 and correlation coefficient r = 0,9908 was obtained. Using this equation the estimated total systematic error of - 0,50% for a level of 150 μg/dL and - 0,75% for a level of 300 μg/dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained:

Samples (µg/dL)	Repetitions	Within-Run Accuracy		Total Accuracy	
		SD (µg/dL)	%CV	SD (µg/dL)	%CV
82,661	80	1,139	1,4	3,388	4,1
163,876	80	1,231	0,8	2,930	1,8
241,319	80	1,381	0,6	4,223	1,8

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

	µg/dL
MEN	65 - 170
WOMEN	50 - 170

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): µmol / L

Iron (µg/dL) x 0,179 = Iron (µmol/L)

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 620 nm (620-630).
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes.

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS)

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotécnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiration date mentioned on the presentation packaging since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotécnica Ltda, phone +55 35 3214 4646 or by email sac@biotecnicaltda.com.br

AUTOMATION

This is automated in most analyzers. The protocols are available at www.biotechnica.ind.br

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