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ENGLISH
INTENDED USE
 Kit intended to determination of iron in the serum. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 15 to 30 ° C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
 The iron present in the sample, at acidic pH, is released from transferrin. Fe³⁺ ions are reduced by hydroxylamine to Fe²⁺. These react with the chromogen ferrozine to form a colored complex, which exhibits a maximum absorption at 560 nm. The intensity of the color is proportional to the concentration of iron in the sample.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practices. The samples to be analyzed should be treated as potentially infectious material.
Preservation: The serum iron is stable for 7 days if stored at a temperature of 20 to 25 °C, 3 weeks if stored at a temperature of 4 to 8 ° C and 1 year if stored at -20 ° C.

PRODUCT HANDLING

R 1	Succinate buffer> 0.1 mol/L; Hydroxylamine hydrochloride> 0.1 mol/L; detergent; reducing agent.	
R 2	Ferrozine> 1.0 mmol /L; preservative.	
STD	Solution of Fe3 + ions in concentration equivalent to 100 µg /dL, preservative, stabilizer. Traceable to reference material NIST -3126a and NIST-928.	

STABILITY IN USE

- The stability of the product (R1, R2 and STD) in use is 24 months, since it is followed recommended storage conditions (15 to 30 ° C).

PRODUCT TREATMENT AND HANDLING

A) REAGENT PREPARATION
 R1 and R2: Ready-to-use reagents.

B) OPERATING INTERVAL
 The operating range of the product is from 6.66 µg/dL to 500.0 µg/dL. For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor .

QUALITY CONTROL
 The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Control of Laboratory Quality it is recommended to use the calibrator serum and the control sera below:

Serum Calibrator - Autocal H	13.002.00
Normal Control Serum – Quantinorm	13.003.00
Serum Pathological Control– Quantialt	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE
 1. Pipette in the assay tubes:

	White	Reactive White (B1)	Standard (P1)	White Sample (A1)	Sample (A2)
Purified water	450 µL	250 µL	----	200 µL	----
Sample	----	----	----	250 µL	250 µL
STD	----	----	250 µL	----	----
R1	800 µL	800 µL	800 µL	800 µL	800 µL
R2	----	200 µL	200 µL	----	200 µL

- Homogenize gently. Incubate for 10 minutes at a temperature of 15 to 30 ° C.
- Measure absorbances: White Reactive (B1), Standard (P1), Sample White (A1) and Sample (A2), resetting the instrument to White. Read absorbances at 560 nm (540-580).

B) CALCULATIONS
Calculation for serum
 Iron (µg/dL) = $\frac{A2 - A1 - B1}{P1 - B1} \times \text{Standard Concentration (µg/dL)}$

Example:
 Standard Concentration = 100 µg/dL
 Absorbance of White Reactive (B1) = 0,03
 Absorbance of standard (P1) = 0,130
 Absorbance of White sample (A1) = 0,02
 Absorbance of sample (A2) = 0,140
 Iron (µg/dL) = $\frac{0,140 - 0,02 - 0,03 \times 100}{0,130 - 0,03} = 90$ µg/dL

With calibration factor:
 Calibration factor = $\frac{\text{Standard Concentration (µg/dL)}}{(P1 - B1)}$
 Iron (µg/dL) = (A2 - A1 - B1) x Calibration factor

Example:
 Calibration factor = $\frac{100}{(0,130 - 0,03)} = 1000$
 Iron (µg/dL) = (0,140 - 0,02 - 0,03) x 1000 = 90 µg/dL

C) INTERPRETATION

In the duodenum and upper jejunum, absorption of iron ingested as Fe²⁺ occurs. The free trivalent form and Fe³⁺ bound to the hemoproteins present in the diet are reduced by vitamin C. About 1 mg of iron is absorbed daily. Once in the mucosal cells, Fe²⁺ ions bind to transporters. Before passing into plasma, these are oxidized by ceruloplasmin to Fe³⁺ and bind to apotransferrin. The apotransferrin-Fe³⁺ complex is called transferrin which carries the iron ions into the blood plasma. The determination of iron (non-heme) is used in the diagnosis and treatment of diseases such as anemia due to iron deficiency, hemochromatosis, chronic kidney disease and in the monitoring of other types of anemias.

INTERFERING AND LIMITATIONS

Hemolysis, jaundice and Lipemia: Hemoglobin > 200 mg/dL / Bilirubin > 40 mg/dL / Triglycerides > 900 mg/dL interferes with the dosage.

Medications: see recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 4.24 µg/dL / Limit of quantification: 6.66 µg/dL

Analytical Specificity: The product specifically determines iron in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 1.002x - 0.165$ and correlation coefficient $r = 0.9999$ was obtained. Using this equation the estimated total systematic error of 0.035% for a level of 100 µg/dL and 0.152% for a level of 350 µg/dL.

Accuracy: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained:

Samples (µg/dL)	Repetitions	Within-Run Accuracy		Total Accuracy	
		SD (µg/dL)	%CV	SD (µg/dL)	%CV
64,780	80	0,711	1,1	1,932	3,0
222,298	80	3,145	1,4	6,092	2,7
418,841	80	1,233	0,3	6,640	1,6

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator, and reagent.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to its use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are key factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

	µg/dL
MEN	65 - 170
WOMEN	50 - 170

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.
Conversion to International System Unit (SI): µmol / L
Iron (µg/dL) x 0.179 = Iron (µmol/L)

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 560 nm (540-580).
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes.

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Safety Data Sheet for this product, available at www.biocentria.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS)

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotécnica reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiration date mentioned on the presentation packaging since they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotécnica Ltda, phone +55 35 3214 4646 or by email sac@biocentrialdta.com.br

AUTOMATION

This is automated in most analyzers.

The protocols are available at www.biocentria.ind.br

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