LIPASE

LIPASE / LIPASA
Ref. 11.009.00

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ENGLISH

INTENDED USE

Kit for the quantitative determination of lipase enzymatic activity in serum and plasma. In vitro diagnostic use.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Lipase catalyzes the cleavage of the 1-2-O-dilauryl-rac-glycerol-3-glutaric-(6'-methylresorufin)-ester substrate in alkaline medium forming 1-2-O-dilauryl-rac-glycerol and glutaric-6'-methylresorufin ester which spontaneously decomposes into glutaric acid and methylresorufin. The velocity of formation of methylresorufin is proportional to the catalytic activity of the lipase in the sample.

1-2-O-dilauryl-rac-glycerol-(6'-methylresorufin)-ester $\xrightarrow{\text{Lipase}}$
1-2-O-dilauryl-rac-glycerol + Glutaric-6'-methylresorufin-ester (unstable)
→ Glutaric Acid + Methylresorufin

SAMPLES – TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum and plasma (heparin).

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Serum and plasma lipase is stable for 7 days at 20-25 °C, 3 weeks at 4-8 °C and 1 year at -20 °C.

PRODUCT DESCRIPTION

R 1	Bicine buffer \geq 50 mmol/L; Colipase \geq 1 mg / L; sodium taurodeoxycholate \geq 5 mmol/L, sodium deoxycholate \geq 1 mmol/L, activator, detergent, preservative. Tartarato buffer \geq 2 mmol/L; sodium taurodeoxycholate \geq 5 mmol/L, Colipase \geq 0,1 mg / L; 1-2-O-dilauryl-rac-glycerol-3-glutaric-(6'-methylresorufin)-ester \geq 2 mmol/L; Stabilizers; activators, detergents, preservative.	X
R 2		X

Traceable to the specific absorbance of the chromogen Methylresorufin using calibrated pipettes and handheld spectrophotometer that provide absolute values.

STABILITY IN USE

- After opening, the product (R1 and R2) in use is stable to the expiration date printed on the label, recommended storage conditions are followed (2 to 8 °C).
- Reagents should remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) REAGENT PREPARATION

R1 and R2: Ready-to-use reagents.

B) OPERATING INTERVAL

The operating range of the product is 3.78 U / L at 300 U / L.

For higher values, dilute the sample with 150 mM NaCl (0.9%), carry out a new dosage and multiply the result obtained by the dilution factor.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and control sera below:

Serum Calibrator - Autocal H	REF	13.002.00
Normal Control Serum- Quantinorm		13.003.00
Serum Pathological Control - Quantialt		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

- 1- Adjust the temperature of the photometer to 37 °C, select the wavelength at 580nm and clear with purified water.
- 2 - Pipette into test tube.

	Volume
R1	800 µL
Purified Water / Sample / CAL	10 µL
R2	200 µL

- 3 - Homogenize and insert in the thermostated cuvette. Set the stopwatch.

- 4- Note the absorbance at 90 seconds (A₁) and 180 seconds (A₂), white, sample and calibrator.

B) CALCULATIONS

Lipase (U / L) = $\frac{\Delta A \text{ Sample} - \Delta A \text{ White}}{\Delta A \text{ Calibrador} - \Delta A \text{ White}} \times \text{Concentration Calibrator (U / L)}$

Where:

$\Delta A \text{ White} = A_2 - A_1$ of purified water.

$\Delta A \text{ Sample} = A_2 - A_1$ of the sample.

$\Delta A \text{ Calibrador} = A_2 - A_1$ of the CAL.

C) INTERPRETATION

Lipase is an enzyme that hydrolyzes esters of glycerol with long chain fatty acids. Most of the enzyme in the serum is derived from the pancreas. Its determination is used to investigate pancreatic disorders, such as acute pancreatitis, chronic pancreatitis, and pancreatic canal obstruction. After an acute pancreatitis event, serum lipase activity increases within 4 to 8 hours, peaking at around 24 hours and decreasing after 8 to 14 days. However, elevation of serum lipase activity is not proportional to the damage caused by the event. Values are also elevated in acute and chronic renal disease and in acute cholecystitis.

INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin> 150 mg / dL / Bilirubin> 40 mg / dL / Triglicérides> 900 mg / dL.

Medications: consult recommended literature reference (Young, 2000).

PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 3.16 U / L / Limit of quantification: 3.78 U / L.

Analytical Specificity: The product specifically determines lipase in the presence of other interfering substances in the sample up to the concentrations reported above.

Accuracy: The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation $y = 0.989x + 2.529$ and correlation coefficient $r = 0.9976$ was obtained. Using this equation the estimated systematic error of 3.958% for a level of 50 U / L and 1.429% for a level of 100 U / L.

Precision: It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, being obtained.

Samples (U/L)	Repetitions	In-run Precision		Total Precision	
		SD (U/L)	%CV	SD (U/L)	%CV
31,480	80	0,540	1,7	0,794	2,5
85,252	80	1,405	1,6	2,540	3,0
147,311	80	2,199	1,5	4,017	2,7

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange reagent bottle caps to avoid cross-contamination.
- Do not use the reagent when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the reagents out of the specified storage conditions when they are not in use.
- Use specific glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The water level of the water bath should be higher than that of the test tubes containing the reactions.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

REFERENCE RANGES

Serum / Plasma	< 60 U/L
----------------	----------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.
Conversion to International System Unit (SI): $\mu\text{kat} / \text{L}$
Lipase (U / L) x 0.0167 = Lipase ($\mu\text{kat} / \text{L}$)

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Spectrophotometer or photometer for reading at 580 nm (550-600).
- Water bath, thermostated at 37 ° C.
- Glass and / or automatic pipettes.
- Clock or Stopwatch.
- Test tubes

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

AUTOMATION

This product is compatible to the most types of biochemical automatic analysers. The applications are available at www.biotechnicaltda.ind.br

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