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### ENGLISH

#### INTENDED USE

Kit intended to determination of total proteins in serum. Use in diagnostics in vitro.

#### STORAGE AND HANDLING

- Store at room temperature.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

#### WORKING PRINCIPLE

The Cu<sup>2+</sup> ions, in alkaline medium, react with the peptide bonds of the proteins forming a color complex that shows maximum absorption at 550 nm. The color intensity is proportional to the protein concentration in the sample.

#### SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION

Sample Type: Serum.

Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practices. The samples to be analyzed should be treated as potentially infectious material.

Preservation: Serum proteins are stable for 3 days if stored at 4 to 8 °C and 7 days if stored at -20 ° C.

#### PRODUCT DESCRIPTION

<b>R 1</b>	Copper sulphate ≥ 5 mmol /L, Sodium potassium tartrate ≥ 20 mmol /L, potassium iodide ≥ 10 mmol/L, Sodium hydroxide ≥ 0.1 mol/L, detergent.	
<b>STD</b>	Phosphate Buffer ≥ 20 mmol / L; Bovine albumin at a concentration equivalent to 5.0 g / dL; Preservative. Traceable to NIST reference material 927d.	

### STABILITY IN USE

- The stability of the product (R1 and STD) in use is 24 months since recommended storage conditions are followed (room temperature).

### TECHNICAL PROCEDURE

#### A) REAGENT PREPARATION

**R1 and STD: Product** is ready to use.

#### B) OPERATING RANGE

The product operating range is from 0,28 g/dL to 12,00 g/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

#### QUALITY CONTROL

The use of controls should be a routine practice in the laboratory. It is suggested to use one control in the reference range or in the decision level and another control with value in another range of clinical significance. For Calibration and Internal Control of Laboratorial Quality it is recommended to use the calibrator serum and the control sera below:

Serum Calibrador - Autocal H	REF	13.002.00
Normal Control Serum - Quantinorm		13.003.00
QuantinormSerum Pathological Control - Quantialt		13.004.00

### TEST PROCEDURE, CALCULATION AND INTERPRETATION

#### A) TEST PROCEDURE

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R1	1.0 mL	1.0 mL	1.0 mL

2. Homogenize and keep for 2 minutes at room temperature.
3. Measure the absorbance of Standard and Sample against the blank at 550 nm. The color is stable for 2 hours.

#### B) CALCULATIONS

Protein (g / dL) =  $\frac{\text{Sample Absorbance} \times \text{STD Concentration (g / dL)}}{\text{STD Absorbance}}$

Example:

STD concentration = 5 g/dL

Sample Absorbance = 0.352

STD Absorbance = 0.250

Protein (g / dL) =  $\frac{0,352 \times 5}{0,250} = 7.04$  g/dL

#### With Calibration Factor:

Calibration Factor =  $\frac{\text{STD Concentration (g / dL)}}{\text{STD Absorbance}}$

Protein (g/dL) = Sample Absorbance x Calibration Factor

Example:

Calibration Factor =  $\frac{5}{0,250} = 20$

Protein (g / dL) = 0,352 x 20 = 7,04 g / dL

#### C) INTERPRETATION

Proteins are predominantly synthesized in the liver, plasma cells, lymph nodes, spleen and bone marrow. In the course of a disease, the total protein concentration and its fractions may deviate significantly from the normal value. Measurement of the total protein can be used in the diagnosis and treatment of a variety of diseases involving the liver, kidney or bone marrow as well as in other metabolic or nutritional disorders. Hypoproteinemia can be caused by diseases and disorders such as: blood loss, nephrotic syndrome, severe burns, salt retention syndrome and Kwarshiorkor (acute protein deficiency). Hyperproteinemia can be observed in cases of severe dehydration and in diseases such as multiple myeloma. Alterations in the relative percentage of plasma proteins can occur without change in the amount of total protein, being caused by variation in the percentage of one of its fractions. The A / G ratio is commonly used as an index of distribution of the albumin and globulin fractions. Changes in this ratio can be observed in liver cirrhosis, glomerulonephritis, nephrotic syndrome, acute hepatitis, lupus erythematosus as well as in certain acute and chronic inflammations.

#### INTERFERING AND LIMITATIONS

Hemolysis, Jaundice and Lipemia: Hemoglobin > 150 mg/dL / Bilirubin > 40 mg/dL / Triglicérides > 288 mg/dL interferes with the sample.

Medications: consult recommended literature reference (Young, 2000).

## PERFORMANCE CHARACTERISTICS

Sensitivity: Detection limit: 0.16 g/dL / Limit of quantification: 0.28 g/dL  
Analytical Specificity: The product specifically determines the protein in the presence of other interfering substances in the sample, up to the concentrations reported above.

**Accuracy:** The method was compared with a similar method by the determination of 40 samples in duplicate. The regression equation  $y = 0.991x + 0.037$  and correlation coefficient  $r = 0.9974$  was obtained. Using this equation the total systematic error estimated from 0.33% for a level of 3.0 g/dL and -0.48% for a level of 9.0 g/dL.

**Accuracy:** It was determined using samples at 03 decision levels, two runs in duplicate per day for 20 days, obtaining:

Samples (g/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD (g/dL)	%CV	SD (g/dL)	%CV
3,6991	80	0,0421	1,1	0,0475	1,3
6,6642	80	0,0348	0,5	0,0450	0,7
10,2555	80	0,0410	0,4	0,0410	0,4

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

## RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

## REFERENCE RANGES

Serum	6,3 - 8,3 g/dL
-------	----------------

These values are for guidance only and it is recommended that each laboratory establishes its own reference range.

Conversion to International System Unit (SI): g/L

Total Protein (g/dL) x 10 = Total Protein (g/L)

## ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer able to read at 550 nm (540-560).
- Glass pipettes and/or automated.
- Clock or chronometer.
- Assay tubes.

## WARNINGS AND PRECAUTIONS FOR DISPOSAL OF PRODUCTS

Disposal, Safety and First Aid information are described in the Material Safety Data Sheet for this product, available at [www.biotechnica.ind.br](http://www.biotechnica.ind.br) or by telephone (35) - 3214-4646.

Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

## QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being released for consumption, all Biotechnique reagents are tested by the Quality Control Department. The quality of the reagents is ensured by the expiration date mentioned on the presentation packaging once they are stored and transported under specified conditions. The data related to the Quality Control of this product (lot printed on the labels of the reagent bottles) or any doubt in the use of this kit, contact the Scientific Advisory of Biotécnica Ltda, by phone +55 35 3214 4646 or by email [sac@Biotechnicaltda.com.br](mailto:sac@Biotechnicaltda.com.br)

## AUTOMATION

This procedure is automated in most analyzers. The protocols are available at [www.biotechnicaltda.ind.br](http://www.biotechnicaltda.ind.br)

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