QUANTINORM
QUANTINORM / QUANTINORM
Ref. 13.003.00

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ENGLISH

INTENDED USE

Kit for the internal quality control of quantitative tests in biochemistry. Use *in vitro* diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Human lyophilized serum used in the control of biochemistry assays.

PRODUCT DESCRIPTION

Human serum plus tissue extracts of human and / or animal origin.

The control shows negative results for the HBs antigen for the HIV and HCV antibody. However, they should be treated as potentially infectious.

The control values are traceable to reference methods and reference materials from the National Institute of Standard and Technology (NIST) and European Reference Materials (ERM).

STABILITY IN USE

- After reconstitution of the control, the analytes have the following stability:
 - Bilirubins and Phosphatase Alkaline: 48 hours when stored at 2 to 8 °C in the original bottle.
 - Other analytes: 5 days when stored at 2 to 8 °C in the original bottle; and 15 days at -20 °C in *eppendorf* flask.
- The control must remain outside the specified temperature only for the time required to perform the tests.

TREATMENT AND HANDLING OF THE PRODUCT

A) CONTROL PREPARATION

Reconstitute with 5.0 mL of purified water. Homogenize gently by inversion, avoiding foaming until complete dissolution of the control.

QUALITY CONTROL

The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For internal control of laboratory quality it is recommended to use the controls below:

QUANTINORM Control	REF	13.003.00
QUANTIALT Control		13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

The control should be used in the same way as patient samples in accordance with the instructions for use of the corresponding biochemical method.

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use PPE according to Good Clinical Laboratory Practice.
- Do not use the calibrator when it has a visual characteristic that is not in accordance with the product's MSDS requirements.
- Avoid leaving the calibrators out of the specified storage conditions when they are not in use.
- Use specific disposable tips for each sample, calibrator and control.
- The laboratory shall establish the chemical, microbiological and particulate requirements for water prior to use for each of its applications and shall define the specifications or types of water that meets them. Once the required purity has been defined, the purification system must be validated and it is important to ensure that the water obtained continues to meet the specifications by means of periodic checks.
- Proper cleaning and drying of the used material are fundamental factors for the stability of the reagents and obtaining correct results.

MATERIALS REQUIRED TO CARRY OUT THE TEST

- Glass and / or automatic pipettes.
- Purified Water

ALERTS AND PRECAUTIONS REGARDING PRODUCT DISPOSAL

- Disposal, Safety and First Aid information are described in the Individual Product Safety Data Sheet (MSDS) for this product, available at www.biotechnica.ind.br or by telephone (35) -3214-4646.
- Dispose of leftover reactions in accordance with Good Clinical Laboratory Practice (BPLC) and Health Service Waste Management Program (PGRSS).

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE

Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this

procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotechnicaltda.com.br

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- CEMBROWSKI, G. Thoughts on quality-control systems: A laboratorians perspective. *Clin. Chem.* v.43, p.886-892, 1997.
- YOUNG, D.S. *Effects of drugs on clinical laboratory tests - vol. 2*, 5 ed. Washington DC: AACCP Press, 2000.
- WESTGARD, J. O. et al. A multi-rule shewhart chart quality control in clinical chemistry. *Clin. Chem.* v.27 p.493-501, 1981. Approved Guideline. *CLSI document H47-A2*. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.

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The values reported for the control serve as a reference. Each laboratory should establish the range of variation in its own analytical system.

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