	Responsável Técnico: Dr. Gilson Sérgio Pizzo CRF MG – 5310 MS 80027310189
		Triglicérides Triglycerides Triglicéridos Ref. 10.010.00

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ENGLISH

INTENDED USE




Kit intended to determination of triglycerides in serum and plasma. Diagnostic use only.

STORAGE AND HANDLING

- Store at 2 to 8 °C.
- Protect from light.
- Stable until the kit expiration date that is printed on the label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE

Triglycerides are hydrolyzed by lipoprotein lipase producing free glycerol, which is transformed into glycerol-3-phosphate by the action of glycerolkinase. The glycerol 3-phosphate is oxidized to dihydroxyacetone and hydrogen peroxide by glycerol-P oxidase. Hydrogen peroxide reacts with 4-aminoantipyrine and 4-chlorophenol under the catalytic action of peroxidase, producing a pink compound (quinonimine) having an absorption maximum at 505 nm.

SAMPLE – PREPARATION AND STABILITY

Sample Type: Serum and plasma (EDTA).
Collection, handling and preparation: Make the collection of the sample in accordance with Good Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.

Preservation: The triglycerides is stable in serum and plasma (EDTA) for 3 days from 4 to 8 °C and 30 days from -20 °C.

PRODUCT DESCRIPTION

	Pipes buffer ≥ 20 mmol/L; 4-chlorophenol ≥ 1 mmol/L; ; 4-aminoantipyrine ≥ 0,1 mmol/L; ATP - adenosine triphosphate ≥ 0,5 mmol/L; Glycerol kinase ≥ 500 U/L; Peroxidase ≥ 1000 U/L; Lipoprotein lipase ≥ 1000 U/L; Glycerol-3-phosphate oxidase ≥ 1000 U/L; activators; detergents, stabilizers and preservative.	
	Glycerol is equivalent to the triglyceride concentration of 200 mg / dL, preservatives.	

Traceable to NIST reference material 1951b.

STABILITY IN USE

- The stability of the product (R1 e STD) in use is 15 months as long as followed by the recommended storage conditions (2 to 8 °C).
- Reagents must remain outside the specified temperature only the time required for testing

TECHNICAL PROCEDURE

A) REAGENT PREPARATION
Reagent 1 (R1)
 Ready-to-use Reagent

B) OPERATING RANGE
 The product operating range is from 9,00 mg/dL to 800 mg/dL. For higher values, dilute the sample with NaCl 150 mM (0,9%), repeat the assay and multiply the result by the dilution factor.

QUALITY CONTROL

Use of controls should be a routine practice in the lab. It is suggested to use a control within the reference range or decision level and another control with value in another range of clinical significance. For the internal quality control of the lab it's indicated the use of the calibrator serum and control sera below:
 Calibrator Serum - Autocal H 13.002.00
 Normal Control Serum – Quantinorm 13.003.00
 Pathological Control Serum – Quantialt 13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION

A) TEST PROCEDURE

1. Pipette in the assay tubes:

	Blank	Standard	Sample
STD	-	10 µL	-
Sample	-	-	10 µL
R 1	1,0 mL	1,0 mL	1,0 mL

2. Homogenize and incubate 10 minutes at 37 °C.
3. Measure the absorbance of Standard and Sample against the blank at 505 nm (490-510 nm). The final reaction is stable 25 minutes.

B) CALCULATIONS

Calculations
 Triglycerides (mg/dL) = $\frac{\text{Sample Absorbance}}{\text{STD Absorbance}} \times \text{STD concentration (mg/dL)}$

Example:
 Sample Absorbance = 0,277
 STD Absorbance = 0,248
 STD concentration = 200 mg/dL
 Triglycerides (mg/dL) = $\frac{0,277}{0,248} \times 200 = 223 \text{ mg/dL}$

With Calibration Factor (CF):
 $CF = \frac{\text{STD concentration}}{\text{Standard Absorbance}}$
 Triglycerides (mg/dL) = Sample Absorbance x CF
Example:
 $CF = \frac{200}{0,248} = 806,5$
 Triglycerides (mg/dL) = 0,277 x 806,5 = 223 mg/dL

C) INTERPRETATION
 Triglycerides are esters of glycerol. Part is synthesized in the liver and another part in the feed is obtained. The triglycerides from the diet are digested in the duodenum and absorbed in the proximal ileum. Through the action of pancreatic and intestinal lipases in the presence of bile acids, they are first hydrolyzed to glycerol, monoglycerides and fatty acids. After absorption, these components again form triglycerides in the intestinal epithelial cells and then are combined with cholesterol and apolipoproteins to form chylomicrons, which are secreted into lymphatic system to then reach the circulation. Triglycerides are the primary metabolic fuel carried by chylomicrons and are distributed to the liver and peripheral cells which are hydrolyzed into fatty acids by lipases.

The determination of triglycerides is used in the diagnosis and treatment of patients with diabetes mellitus, nephrosis, liver obstruction, disorders in the metabolism of lipids and several other endocrine diseases.

INTERFERING AND LIMITATIONS
Anticoagulants: Citrate, fluoride, heparin and oxalate may interfere in the reaction.
Hemolyzed, jaundice and Lipemic Sera: Hemoglobin > 500 mg/dL, Bilirubin > 35 mg/dL.
Drugs: see recommended reference (Young, 2000).

PERFORMANCE CHARACTERISTICS
Sensitivity: Detection limit: 1,38 mg/dL / quantification limit: 9,00 mg/dL.
Analytical Specificity: The product determines triglycerides specifically in the presence of other interfering substances in the sample, up to the concentrations stated above.
Accuracy: The method was compared with a similar method for the determination of 40 samples in duplicate. The regression equation obtained was $y = 0,995x - 0,950$ and the correlation coefficient $r = 0,9997$. Using this equation the total systematic error estimated is -1,02% to a level of 180 mg/dL and -0,69% to a level of 500 mg/dL.
Precision: Determined using samples in 03 levels of decision, with two runs in duplicate per day during 20 days, being obtained.

Samples (mg/dL)	Repetitions	Within-Run Precision		Total Precision	
		SD	%CV	SD	%CV
104,686	80	1,173	1,1	1,176	1,1
227,773	80	3,212	1,4	3,770	1,7
626,552	80	2,090	0,3	2,228	0,4

% CV: Coefficient of variation expressed as a percentage; SD: Standard Deviation

RESIDUAL RISK, WARNINGS AND PRECAUTIONS

- Use Protective Equipment in accordance with Good Laboratory Practice.
- Do not mix reagents from different lots.
- Do not exchange the caps from different reagents in order to avoid cross contamination.
- Do not use the reagent if it displays any sign in disagreement with the one specified in the product MSDS.
- Avoid leaving reagents outside the specified storage conditions, when they are not in use.
- Use glass pipettes and disposable tips for each sample, control, standard / calibrator and reagent.
- The level of the water bath must be greater than the test tubes containing the reaction.
- The laboratory should establish chemical, microbiological and particulates requirements for the water prior to its use for each of its applications and must define the specifications or types of water that meets them. Once the required purity has been set, the purification system must be validated and it is important to ensure that the resulting water continues to meet specifications by means of periodic checks.
- The cleaning and drying of the lab material are fundamental factors for the stability of the reagents and reaching correct results.

REFERENCE RANGES
 In conformity to the Brazilian Directives for Dislipidemias and Prevention of Atherosclerosis the reference value for Triglycerides may be different according to the age group as shown in the table below:

Age	Value (mg/dL)	Interpretation
from 2 to 19 years old	< 100	Desirable
	100 – 129	Limitrophe
	> 130	High
More than 20 years old	< 150	Desirable
	150 – 200	Limitrophe
	200 - 499	High
	≥ 500	Very High

These values are intended for orientation only. It's recommended that each lab establishes its own reference range.

Conversion for International System of Units (SI): mmol/L
 Triglycerides (mg/dL) X 0.0113 = Triglycerides (mmol/L)

ADDITIONAL EQUIPMENT

- Spectrophotometer or photometer with a thermostated cuvette able to read at 505 nm (490 – 510 nm).
- Thermostatic water bath at 37°C.
- Glass Pipettes and/or automatic.
- Clock or chronometer.
- Assay tubes

WARNINGS AND PRECAUTIONS

- Discard the reactions surplus according to Good Laboratory Practice in a proper place for potentially infective material.
- The information for Disposing, Security and First Aid are described in the Manual Safety Data Sheet (MSDS) of this product available at www.biotecnica.ind.br or calling for 55-35-3214-4646

QUALITY ASSURANCE / CUSTOMER TECHNICAL SERVICE
 Before being approved for use BioTécnica reagents are tested in the Quality Control Department. The quality of the reagents is assured up to the expiring date stated in the label of the external packaging, since it is stored and transported in the specified conditions. The quality control data concerning this product (batch printed on the labels of reagent bottles) or any technical doubt on handling this product or this procedure, contact us calling +55 35 3214 4646, your local distributor or sending an e-mail: sac@biotecnicaalda.com.br

AUTOMATION
 This product is compatible to the most types of biochemical automatic analyzers. The applications are available at www.biotecnicaalda.ind.br

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