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ENGLISH
INTENDED USE
Kit for the determination of urea in serum, plasma and urine. Use in in vitro diagnostic.

STORAGE AND HANDLING

- Store at 2 to 8 ° C.
- Protect from light.
- The validity of the kit is printed on the package label.
- Do not use reagents whose date has expired.

WORKING PRINCIPLE
The urea of the sample is hydrolyzed by urease with production of carbonic gas and ammonium ions. It reacts with salicylate and sodium hypochlorite in alkaline medium, in the presence of nitropirussiato, producing indophenol, green in color, which can be measured in a spectrophotometer at 580 nm. The intensity of the color is proportional to the concentration of urea in the sample.

SAMPLES: TYPE, COLLECTION, HANDLING, PREPARATION AND PRESERVATION
Sample Type: Serum, plasma (EDTA and heparin) and urine.
Collection, handling and preparation: Perform the sample collection according to Good Clinical Laboratory Practice. The samples to be analyzed should be treated as potentially infectious material.
Preservation: Urea in serum and plasma (EDTA and heparin) is stable for 7 days if stored at 20 to 25 ° C, 7 days if stored at 4 to 8 ° C and 12 months at -20 ° C. Urea in the urine is stable for 2 days if stored at 20 to 25 ° C, 7 days if stored at 4 to 8 ° C and 28 days if stored at -20 ° C. Collect 24-hour urine in a vial containing 2 mL of HCl 50%. Centrifuge for 10 minutes at 3000 rpm. Dilute an aliquot of the urine centrifuged at the ratio of 1:50 with purified water and perform the assay. Multiply the result obtained by 50. This dilution should result in a measurement within the operating range. For higher values perform new dilution by changing the ratio.

PRODUCT DESCRIPTION

R 1	Phosphate buffer ≥ 3 mmol/L; Sodium salicylate ≥ 10 mmol/L; Sodium nitroprusside ≥ 1 mmol/L; stabilizers; preservative.	
R 2	Sodium hypochlorite ≥ 0,01% v/v; Sodium hydroxide ≥ 100 mmol/L.	

R 3	Phosphate buffer ≥ 10 mmol/L; Glycerol ≥ 20% v/v; Urease ≥ 50.000 U/L; stabilizers; preservative.	
STD	Phosphate buffer ≥ 2 mmol/L; Urea at a concentration equivalent to 70 mg/dL; preservative. Traceable to NIST reference material 912a.	

STABILITY IN USE

- After opening, the product (R1, R2, R3 and STD) in use is stable up to the expiration date printed on the label since recommended storage conditions (2 to 8 ° C) are followed.
- The stability of the working reagent is 21 days, since recommended preparation and storage conditions are followed (2 to 8 ° C).
- Reagents must remain outside the specified temperature only the time required for testing.

TREATMENT AND HANDLING OF THE PRODUCT
A) REAGENT PREPARATION
Working Reagent (R1 and R3):
Mix in proportion: 25 parts of R1 + 1 part of R3. Homogenize gently.
R2 and STD:
Ready for use.

B) OPERATING INTERVAL
Serum and Plasma: The operating range of the product is from 2,206 mg/dL to 200 mg/dL. For higher values, dilute the sample with 150 mM NaCl (0,9%), carry out a new dosage and multiply the result obtained by the dilution factor.
Urine: The operating range of the product is from 2,303 g/L to 100,00 g/L corrected by the dilution factor.

QUALITY CONTROL
The use of controls should be routine practice in the laboratory. It is suggested to use one control in the reference range or at the decision level and another control with value in another range of clinical significance. For Calibration and Internal Quality Control Laboratory it is recommended to use the calibrator serum and the control sera below:

Calibrator Serum - Autocal H	13.002.00
Normal Control Serum - Quantinorm	13.003.00
Pathological Control Serum - Quantialt	13.004.00

TEST PROCEDURE, CALCULATION AND INTERPRETATION
A) TEST PROCEDURE
1. Pipette in the assay tubes:

STD	White	Standard	Sample
---	---	10 µL	---
Sample	---	---	10 µL
Purified Water	10 µL	---	---
Working Reagent	1,0 mL	1,0 mL	1,0 mL

2. Homogenize well and incubate the tubes for 5 minutes at 37 ° C.

3. Add:

R2	1,0 mL	1,0 mL	1,0 mL
-----------	--------	--------	--------

4. Homogenize well and incubate the tubes for 5 minutes at 37 ° C.

5. Measure Standard Absorbance (STD) and Sample vs. Blank at 580 nm. Color is stable for 15 minutes.

B) CALCULATIONS
Serum / Plasma:
Urea (mg / dL) = $\frac{\text{Sample Absorbance}}{\text{Absorbance of STD}} \times \text{STD Concentration (mg / dL)}$

Exemplo:
STD concentration = 70 mg / dL
Sample Absorbance = 0.209
Absorbance of STD = 0.378
Urea (mg / dL) = $\frac{0,209 \times 70}{0,378} = 38,7$ mg / dL

With Calibration Factor:
Calibration Factor = $\frac{\text{STD Concentration (mg / dL)}}{\text{Absorbance of STD}}$


Urea (mg / dL) = Sample Absorbance x Calibration Factor

Exemplo:
Calibration Factor = $\frac{70}{0,378} = 185$
Urea (mg / dL) = 0.209 x 185 = 38.7 mg / dL

Urine:
Determine the concentration of urea in the urine sample using the same calculation procedure for serum. Multiply the value obtained by 50 to obtain the result in mg / dL. This should be used to calculate the value of urea in g / 24h, according to the equation below:
Urea (g / 24 hours) = $\frac{\text{mg / dL} \times \text{urinary volume (in L)}}{100}$

C) INTERPRETATION
Urea is the major nitrogen metabolite of protein catabolism. More than 90% of the urea is excreted through the kidneys, consequently, kidney disease is associated with accumulation of urea in the blood. The dosage of blood and urine urea can provide useful clinical information. It may be elevated because of different factors such as high protein diet, high protein catabolism, blood protein reabsorption after a gastrointestinal bleeding, treatment with cortisol and its analogs, dehydration and reduced kidney perfusion. It is also elevated in obstructive post-renal conditions such as malignant tumors, nephrolithiasis and prostatism.

INTERFERING AND LIMITATIONS
Hemolysis, Jaundice and Lipemia:

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